

**FLOWSAFE II** <sup>®</sup>  
Disposable BiLevel CPAP System

## BiLevel & CPAP Disposable System in One

Flow-Safe II+<sup>®</sup> is the only available Disposable BiLevel CPAP System! The lightweight disposable feature allows for easy CPAP or BiLevel CPAP therapy set up and delivery during transport. All devices include an integrated manometer for verifying CPAP, IPAP and EPAP pressures. The Disposable Flow-Safe II+ System is ideal for situations where backup BiLevel CPAP equipment is scarce or unavailable.



- First Disposable BiLevel Device with CPAP or BiLevel Mode switch allows easy transition from CPAP to BiLevel Mode
- No Extra Parts – Flow Determines CPAP or BiLevel IPAP Pressure, Easy EPAP Dial adjusts expiratory pressure
- Verify CPAP & BiLevel Pressures with Integrated Manometer (available on all models)
- Easy to transport
- No Capital Equipment – Ideal for Acute Care and Pre-hospital Settings



**FLOWSAFE II** +<sup>®</sup>  
Disposable BiLevel CPAP System

**How the System Works**

By connecting the Flow-Safe II<sup>®</sup> BiLevel CPAP System to a flowmeter/regulator start with the green knob in the CPAP position. By adjusting flow rates up and down, various CPAP pressures can be achieved and verified with the manometer.



**Blue knob adjusts EPAP Pressure**

**Green knob selects CPAP or BiLevel pressure mode**

**CONNECT TO FLOW SOURCE ONLY**



Flow (LPM)	CPAP MODE (cm H <sub>2</sub> O)
6	2.0 - 3.0
10	6.0 - 7.0
12	8.0 - 9.0
15	11.0 - 12.0

**CAUTION:** CPAP pressure will decrease when BiLevel is activated & increase when BiLevel is deactivated. Verify CPAP pressure with manometer & adjust flowmeter as needed.

**CONNECT TO FLOW SOURCE ONLY**



Flow (LPM)	BiLevel MODE (cm H <sub>2</sub> O)
14	8 - 9 IPAP
15	9 - 10 IPAP
16	11 - 12 IPAP
17 (MAX)	12 - 13 IPAP

**CAUTION:** CPAP pressure will decrease when BiLevel is activated & increase when BiLevel is deactivated. Verify CPAP pressure with manometer & adjust flowmeter as needed.

For BiLevel therapy, simply switch the green knob to the BiLevel position. The CPAP pressure now becomes the Inspiratory Positive Airway Pressure (IPAP). Adjust the blue EPAP knob to obtain desired Expiratory Positive Airway Pressure.

**Flow-Safe II<sup>®</sup> Disposable CPAP Therapy System Compared to Non-invasive Mechanical Ventilation**



Is the flow-safe disposable continuous positive airway pressure (CPAP) system as effective as non-invasive mechanical ventilation (NIMV) in the treatment of acute cardiogenic pulmonary oedema?

Ilhan UZ, MD<sup>a,\*</sup>, Güçlü Selahattin KIYAN, MD<sup>a</sup>, Enver ÖZÇETE, MD<sup>a</sup>, Sercan YALÇINLI, MD<sup>a</sup>, Mehmet Birkan KORGAN, MD<sup>a</sup>, Yusuf Ali ALTUNCI, MD<sup>a</sup>, Murat ERSEL, MD<sup>a</sup>, Funda Karbek AKARCA, MD<sup>a</sup>, Oğuz YAVUZGİL, MD<sup>b</sup>

<sup>a</sup> Department of Emergency Medicine, Ege University Faculty of Medicine, Izmir, Turkey  
<sup>b</sup> Department of Cardiology, Ege University Faculty of Medicine, Izmir, Turkey

According to this recent American Journal of Emergency Medicine, For acute cardiogenic pulmonary oedema, (ACPO) patients, the Flow-Safe Disposable CPAP device is "as effective as NIMV, in improving blood pressure, pulse, respiratory rate and blood gas."

Flow-Safe II, Flow-Safe II EZ<sup>®</sup> and Flow-Safe II+ can be solutions for situations where non-invasive mechanical ventilation equipment is scarce or unavailable.

ORDERING INFORMATION		
Part #	Description	Packaged
#10-57400	Flow-Safe II+ BiLevel CPAP w/ Large Adult Full Face Mask w/ Straight Swivel Port & Headstrap	5/Box
#10-57401	Flow-Safe II+ BiLevel CPAP w/ Medium/Small Adult Full Face Mask w/ Straight Swivel Port & Headstrap	5/Box
#10-57402	Flow-Safe II+ BiLevel CPAP w/ Small/Child Full Face Mask w/ Straight Swivel Port & Headstrap	5/Box
#10-57403	Flow-Safe II+ BiLevel CPAP w/o Mask	5/Box