<b>MercuryMedical</b> 	<b>(en)</b>
11300 - 49th Street North Clearwater, Florida 33762-4807 USA www.mercurymed.com	

# FLOWSAFE II®

DISPOSABLE CPAP SYSTEM

#### DESCRIPTION:

The Mercury Flow-Safe Continuous Positive Airway Pressure (CPAP) device is a respiratory aid intended for use with a facemask and gas supplying device to elevate pressure in the patient's lungs.

#### CONNECTIONS:

- Standard oxygen tubing nipple
- Patient connection ISO 5356-1 - 15mm taper female and 22mm taper male

#### INDICATIONS FOR USE:

The Mercury Flow-Safe CPAP device is intended to provide CPAP to spontaneously breathing patients in the hospital and pre-hospital environment.

#### WARNINGS:

- Do not allow smoking or use unit near sparking equipment, open flame, oil or other flammable chemicals.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

#### CAUTIONS:

- Federal (USA) law restricts this device to sale by or on the order of a physician.
- Do not clean, soak, rinse or sterilize.
- Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended.
- In the event of undesirable flow rate from oxygen source, simply remove the device and place on supplemental oxygen per protocol.
- Use of the Flow-Safe with non-back pressure compensated flow devices may affect input gas liter flow. Always verify delivered CPAP pressure on manometer.

#### CONTRAINDICATIONS:

- Respiratory Arrest
- Unconscious
- Cardiogenic Shock
- Pneumothorax
- Facial Anomalies
- Facial Trauma

#### DIRECTIONS FOR USE:

- Connect O<sub>2</sub> tubing nipple to gas source.
- Secure the face mask snugly to patient's face using head harness.
- Slowly increase gas flow to 6 or 8 LPM. Check mask fit to patient and device connections for leaks.
- Adjust the flowmeter until desired pressure is obtained. Flow of 12 - 14 LPM is required to reach CPAP pressure of 8.5 - 10 cm H<sub>2</sub>O.
- Do not exceed 30 LPM.
- Patient SaO<sub>2</sub> should be monitored using a pulse oximeter.

#### MEASURING PRESSURE:

- Pressure relief limits maximum CPAP pressure to 25 cm H<sub>2</sub>O @ 25 LPM.
- Do not exceed pressure limit of manometer (25 cm H<sub>2</sub>O).
- Manometer accuracy ± 3 cm H<sub>2</sub>O up to 15 cm H<sub>2</sub>O and ± 5 cm H<sub>2</sub>O over 15 cm H<sub>2</sub>O.

						
2797 <b>EMERGO EUROPE</b> Prinsessegracht 20 2514 AP The Hague The Netherlands	Do Not Reuse	Non Sterile	DINP	Not Made with Natural Rubber Latex	By Prescription Only	MR Conditional

© Registered U.S. Trademark of Mercury Enterprises, Inc.  
Patent # US 8,522,618 B1, US 9,370,635 B2, Other Patents Pending

Made in Malaysia	1/2020	#83-900-0390	Rev. 6
------------------	--------	--------------	--------

<b>MercuryMedical</b> 	<b>(de)</b>
11300 - 49th Street North Clearwater, Florida USA 33762-4807 USA www.mercurymed.com	

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