FLWSAFE II

Disposable BiLevel CPAP System

DESCRIPTION:
The Mercury Medical Flow-Safe II™ device is a respiratory aid intended for use with a facemask or advanced airway and gas supplying device to elevate pressure in the patient's lungs. The level of airway pressure delivered to the respiratory system may be set to increase during inspiration and decrease during expiration.

- CONNECTIONS:**
- Standard oxygen tubing nipple
 - Patient connection ISO 5356-1 – 15mm taper female and 22mm taper male

INDICATIONS FOR USE:
The Mercury Medical Flow-Safe II™ device is intended to provide CPAP or BiLevel CPAP pressures to spontaneously breathing patients in the hospital, surgery center and ambulatory setting.

WARNINGS:

- Do not allow smoking or use unit near sparking equipment, open flame, oil or other flammable chemicals.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

CAUTIONS:
Federal (USA) law restricts this device to sale by or on the order of a physician. Do not clean, soak, rinse or sterilize. Do not attempt to disinfect any part of this device. Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended. In the event of undesirable flow rate from oxygen source, simply remove the device and place on supplemental oxygen per protocol or guidelines for other accepted procedures.

Use of Flow-Safe II™ device with non-back-pressure compensated flow devices may affect input gas flow.
Do not use with a vented mask or mask with anti-asphyxia valve. Recommend use with Mercury Medical **straight port full face deluxe mask**.
Device is not intended for long term use.
If device becomes damaged, sealed or occluded, discard the device.
Never wait to begin CPAP or BiLevel CPAP therapy. If a Flow-Safe II™ system is not immediately available, cannot be used effectively or in the event of an undesirable performance, follow your department protocol or guidelines for other accepted procedures.
To be used only by clinicians who have been trained in BiLevel/CPAP therapy. Adjusting input flow rate will affect CPAP and IPAP/EPAP pressures. ALWAYS verify CPAP and IPAP/EPAP pressures with manometer.
This device is assembled with PVC components which include low levels of a phthalate as indicated. Consideration should be used when treating children or pregnant or nursing women.
Listed values are nominal only. Actual results may vary on device and accuracy of the flowmeter.
This device is not intended to be used as a transport ventilator, and is not intended for life support.
Patient SaO₂ should be monitored using a pulse oximeter.
Patient ECG/Co₂ should be monitored using capnography.
Use in CPAP mode if used in conjunction with a rebuilder.
Recommended therapy should be started in the default CPAP mode.

CONTRAINDICATIONS:

- Respiratory Arrest
- Unconscious
- Cardogenic Shock
- Pneumothorax
- Facial Anomalies
- Facial Trauma
- Airway Obstruction


DIRECTIONS FOR PROVIDING CPAP:
Ensure device is set in CPAP mode. Green switch is set to CPAP (See Figure 1).
1. Connect O₂ tubing nipple to gas source. Turn on gas source.
2. Secure the face mask snugly to patient's face using head harness.
3. Slowly increase gas flow to 6 - 8 LPM. Check mask fit to patient and device connections for leaks.
4. Adjust the flowmeter until desired pressure is obtained. Flow of 12 - 14 LPM is required to reach CPAP pressure of 8.5 - 10 cm H₂O.

DIRECTIONS FOR PROVIDING BiLEVEL:
Ensure device is set in BiLevel mode by rotating the Green switch to the BiLevel setting (See Figure 2).
1. Connect O₂ tubing nipple to gas source. Turn on gas source.
2. Secure the face mask snugly to patient's face using head harness.
3. Slowly increase gas flow to reach approximately 8 cm H₂O IPAP. This is the minimum IPAP pressure required for device to function properly in the BiLevel mode. Check mask fit to patient and device connections for leaks.
4. Adjust the flowmeter until desired IPAP pressure is obtained. Flow of 17 LPM is required to reach maximum IPAP pressure of 12 - 13 cm H₂O.
5. To decrease EPAP pressure rotate EPAP knob counterclockwise. To increase EPAP pressure rotate EPAP knob clockwise. (See Figure 3)
6. Effective mask seal is required for device to shift into EPAP mode from IPAP mode.

MEASURING PRESSURE:

- Internal pressure limits maximum CPAP pressure to 25 cm H₂O @ 25 LPM.
- Do not exceed pressure limit of manometer (25 cm H₂O).
- Manometer accuracy +/-3 cm H₂O up to 15 cm H₂O and +/- 5 cm H₂O over 15 cm H₂O.

MR Conditional: Static magnetic field of 3.0 Tesla or less
Spatial Gradient of 3000 gauss/cm or less

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