(For use with FLOWSAFE FLOWSAFE II FLOWSAFE II EZ FLOWSAFE II-⁺)

INSTRUCTIONS FOR USE:
1. Connect each leg of WYE tubing to gas sources, providing flow ranges.
2. Connect Flow-Safe to WYE.
3. Adjust each flow device for desired flow range. **NOTE:** Flows from each device provide combined flow on outlet.

WARNINGS:
1. Output of WYE flows are combined.
2. If using blended gases (i.e. air/O₂), see chart below:

Oxygen Flowmeter Setting LPM	Air Flowmeter Setting LPM	Oxygen Output Percent
1	29	24
2	28	26
3	27	29
4	26	31
5	25	34
6	24	37
7	23	39
8	22	42
9	21	45
10	20	47
11	19	50
12	18	53
13	17	55
14	16	58
15	15	60

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3. Always monitor output flows with oxygen analyzer, to confirm oxygen concentration.
4. To prevent flow restriction or back pressure, do not use additional connectors or swivels.
5. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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NOTE: The listed output is nominal value only, actual output may vary depending upon patient tidal volume, breath rate and the existence of mask leaks. Faster breathing rates lower oxygen concentration while slower rates result in higher oxygen concentration. Large tidal volumes lower oxygen concentration, small tidal volumes result in higher oxygen concentration.

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(For bruk med FLOWSAFE FLOWSAFE II FLOWSAFE II EZ FLOWSAFE II EZ FLOWSAFE II EZ)

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