Manometer

(en)

Single-Patient-Use

Disposable Airway Pressure Manometer

INDICATIONS FOR USE:


The intended use of the Disposable Manometer is to provide a method to monitor airway pressure.

ACCESSORIES:


Mercury Medical In-Line Pressure Manometer Adapter (left or right)
 Mercury Medical Pressure Manometer Stand and Tubing

WARNINGS:

- Do not connect to pressure sources exceeding the limits of the manometer.
- The use of the in-line adapters will increase airway dead space and their use is not recommended if the manometer can be connected directly to a built in manometer port.

- When using the in-line adapter on resuscitators without built in manometer ports, the Mercury Medical Pressure Manometer must be placed between the patient valve and the patient.

- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

CAUTIONS:

- Single use only. Do not sterilize or re-use to reduce the risk of possible transmission of infection.
- During use, the patient's condition must be continuously monitored.
- Federal (USA) law restricts this device to sale by or on the order of a physician.
- Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended.

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