

Adapters/Fittings/Chemetron (NCG) Style



- #69-20100 Chemetron oxygen adapter x 1/4" hose barb
- #69-20101 Chemetron nitrous oxide adapter x 1/4" hose barb
- #69-20102 Chemetron medical air adapter x 1/4" hose barb
- #69-20103 Chemetron vacuum adapter x 1/4" hose barb
- #69-20364 Chemetron vacuum adapter x 5/16" hose barb
- #69-20363 Chemetron WAGD adapter x 1/4" hose barb
- #69-21510 Chemetron WAGD adapter x 5/16" hose barb
- #69-20366 Chemetron carbon dioxide adapter x 1/4" hose barb (CO₂>7%)
- #69-20367 Chemetron oxygen/carbon dioxide adapter x 1/4" hose barb (95% O₂/5% CO₂) (CO₂<7%)



- #69-20104 Chemetron oxygen adapter x 1/8-27 male
- #69-20105 Chemetron nitrous oxide adapter x 1/8-27 male
- #69-20106 Chemetron medical air adapter x 1/8-27 male
- #69-20107 Chemetron vacuum adapter x 1/8 - 27 male
- #69-20368 Chemetron WAGD adapter x 1/8 - 27 male
- #69-21253 Chemetron carbon dioxide adapter x 1/8 - 27 male
- #69-21272 Chemetron oxygen/carbon dioxide adapter x 1/8-27 male (95% O₂/5% CO₂)



- #69-21290 Chemetron WAGD adapter x D.I.S.S. WAGD female
- #69-20108 Chemetron oxygen adapter x D.I.S.S. oxygen female
- #69-20109 Chemetron nitrous oxide adapter x D.I.S.S. nitrous oxide female
- #69-20110 Chemetron medical air adapter x D.I.S.S. medical air female
- #69-20111 Chemetron vacuum adapter x D.I.S.S. vacuum female
- #69-21291 Chemetron carbon dioxide adapter x D.I.S.S. carbon dioxide female
- #69-21292 Chemetron oxygen/carbon dioxide adapter x D.I.S.S. female



- #69-21288 Chemetron carbon dioxide adapter x D.I.S.S. carbon dioxide male
- #69-21289 Chemetron oxygen/carbon dioxide adapter x D.I.S.S. male (95% O₂/5% CO₂)
- #69-20112 Chemetron oxygen adapter x D.I.S.S. oxygen male
- #69-20113 Chemetron nitrous oxide adapter x D.I.S.S. nitrous oxide male
- #69-20114 Chemetron medical air adapter x D.I.S.S. medical air male
- #69-20115 Chemetron vacuum adapter x D.I.S.S. vacuum male
- #69-21287 Chemetron WAGD adapter x D.I.S.S. WAGD male



- #69-21484 Chemetron vacuum adapter - Chemetron inlet to double vacuum D.I.S.S. Female outlets, center mounted
- #69-20263 Chemetron oxygen adapter - Chemetron inlet to double oxygen D.I.S.S. female outlets, center mounted
- #69-20384 Chemetron medical air adapter - Chemetron inlet to double medical air D.I.S.S. female outlets, center mounted



- #69-21485 Chemetron medical air adapter - Chemetron inlet to double medical air D.I.S.S. female outlets, offset left
- #69-20264 Chemetron oxygen adapter - Chemetron inlet to double oxygen D.I.S.S. female outlets, offset left
- #69-21486 Chemetron vacuum adapter - Chemetron inlet to double vacuum D.I.S.S. female outlets, offset left



- #69-20265 Chemetron oxygen adapter - Chemetron inlet to double oxygen D.I.S.S. female outlets, offset right
- #69-21487 Chemetron medical air adapter- Chemetron inlet to double medical air D.I.S.S. female outlets, offset right
- #69-21488 Chemetron vacuum Chemetron inlet to double vacuum D.I.S.S. female outlets, offset right

Adapters/Fittings/chemetron (NCG) Style



- #69-20116 Chemetron oxygen coupler x 1/4" hose barb
- #69-20117 Chemetron nitrous oxide coupler x 1/4" hose barb
- #69-20118 Chemetron medical air coupler x 1/4" hose barb
- #69-20119 Chemetron vacuum coupler x 1/4" hose barb
- #69-21511 Chemetron vacuum coupler x 5/16" hose barb
- #69-20369 Chemetron WAGD coupler x 1/4" hose barb
- #69-20371 Chemetron carbon dioxide coupler x 1/4" hose barb (CO₂ > 7%)
- #69-20372 Chemetron carbon dioxide coupler x 1/4" hose barb (CO₂ < 7%)



- #69-20386 Chemetron carbon dioxide coupler x 1/8-27 male (CO₂<7%)
- #69-20387 Chemetron carbon dioxide coupler x 1/4-18 male (CO₂ > 7%)
- #69-20121 Chemetron oxygen coupler x 1/8-27 male
- #69-20124 Chemetron oxygen coupler x 1/4-18 male
- #69-20122 Chemetron vacuum coupler x 1/8-27 male
- #69-20125 Chemetron vacuum coupler x 1/4-18 male
- #69-20123 Chemetron medical air coupler x 1/8-27 male
- #69-20127 Chemetron medical air coupler x 1/4-18 male
- #69-20126 Chemetron nitrous oxide coupler x 1/4-18 male
- #69-20374 Chemetron nitrous oxide coupler x 1/8-27 male
- #69-20373 Chemetron WAGD coupler x 1/4-18 male
- #69-20375 Chemetron WAGD coupler x 1/8-27 male



- #69-20120 Chemetron oxygen coupler x D.I.S.S. oxygen male
- #69-20376 Chemetron nitrous oxide coupler x D.I.S.S. nitrous oxide male
- #69-20377 Chemetron medical air coupler x D.I.S.S. medical air male
- #69-20378 Chemetron vacuum coupler x D.I.S.S. vacuum male
- #69-20379 Chemetron WAGD coupler x D.I.S.S. male
- #69-20381 Chemetron carbon dioxide coupler x D.I.S.S. carbon dioxide male



- #69-21293 Chemetron oxygen coupler x 1/8-27 female
- #69-21294 Chemetron nitrous oxide coupler x 1/8-27 female
- #69-21295 Chemetron medical air coupler x 1/8-27 female
- #69-21296 Chemetron vacuum coupler x 1/8-27 female

Adapters/Fittings/Chemetron Style



- #69-21273 Chemetron oxygen adapter x 1/8-27 female
- #69-21274 Chemetron nitrous oxide adapter x 1/8 - 27 female
- #69-21275 Chemetron medical air adapter x 1/8 - 27 female
- #69-21276 Chemetron vacuum adapter x 1/8-27 female
- #69-21277 Chemetron WAGD adapter x 1/8-27 female
- #69-21278 Chemetron carbon dioxide adapter x 1/8-27 female
- #69-21279 Chemetron 95% O₂/5% CO₂ adapter x 1/8-27 female



- #69-21280 Chemetron oxygen adapter x 1/4-18 male
- #69-21281 Chemetron nitrous oxide adapter x 1/4-18 male
- #69-21282 Chemetron medical air adapter x 1/4-18 male
- #69-21283 Chemetron vacuum adapter x 1/4-18 male
- #69-21284 Chemetron WAGD adapter x 1/4-18 male
- #69-21285 Chemetron carbon dioxide adapter x 1/4-18 male
- #69-21286 Chemetron 95% O₂/5% CO₂ adapter x 1/4-18 male



- #69-21809 Chemetron oxygen coupler x D.I.S.S. oxygen female
- #69-21850 Chemetron nitrous oxide coupler x D.I.S.S. nitrous oxide female
- #69-21848 Chemetron medical air coupler x D.I.S.S. medical air female
- #69-21849 Chemetron vacuum coupler x D.I.S.S. vacuum female
- #69-21859 Chemetron WAGD coupler x D.I.S.S. WAGD female
- #69-21880 Chemetron carbon dioxide coupler x D.I.S.S. carbon dioxide female



- #69-21254 Chemetron oxygen duplex - wye style
- #69-21265 Chemetron vacuum duplex - wye style
- #69-21266 Chemetron medical air duplex - wye style
- #69-21271 Chemetron nitrous oxide duplex - wye style



- #69-20370 Replacement, rectangular striker pin - for medical air, vacuum* and WAGD
- #69-20385 Replacement, round striker pin - for oxygen and nitrous oxide
- #32-11-90-0002 Replacement, striker screw (not shown)

NOTE: Pins and screw fit original Chemetron adapters only.
*Post 1966 outlets

Adapters/Fittings



- #69-20189* Chemetron oxygen adapter - center mounted
- #69-20190* Chemetron medical air adapter - center mounted
- #69-20191* Chemetron vacuum adapter - center mounted
- #69-21842* Chemetron nitrous oxide adapter - center mounted



- #69-20192* Chemetron oxygen adapter - offset left mounted
- #69-20262* Chemetron vacuum adapter - offset left mounted
- #69-20382* Chemetron medical air adapter - offset left mounted



- #69-20402* Chemetron oxygen adapter - offset right mounted
- #69-20403* Chemetron vacuum adapter - offset right mounted
- #69-20383* Chemetron medical air adapter - offset right mounted

* Couplers spaced 4-1/2" on-center measurement



- #69-20266 Chemetron vacuum adapter, bottle holder type. Hose barb outlet with on/off switch.



- #69-20267 Chemetron vacuum adapter, bottle holder type. Two hose barb outlets with on/off switch.



- #69-20389 Ohmeda vacuum adapter, bottle holder type. Two hose barb outlets with on/off switch.



- #69-20268 D.I.S.S. vacuum adapter, Chemetron bottle holder type. D.I.S.S. female inlet, hose barb outlet with on/off switch.



- #69-20269 D.I.S.S. vacuum adapter, Chemetron bottle holder type. D.I.S.S. female inlet, two hose barb outlets with on/off switch.

Adapters/Fittings/ohmeda (Ohio) Style



- #69-20138 Ohmeda type oxygen adapter x 1/4" hose barb
- #69-20139 Ohmeda type nitrous oxide adapter x 1/4" hose barb
- #69-20140 Ohmeda type medical air adapter x 1/4" hose barb
- #69-20141 Ohmeda type vacuum adapter x 1/4" hose barb
- #69-21512 Ohmeda type vacuum adapter x 5/16" hose barb
- #69-20501 Ohmeda type WAGD adapter x 1/4" hose barb
- #69-21513 Ohmeda type WAGD adapter x 5/16" hose barb
- #69-20760 Ohmeda type carbon dioxide adapter x 1/4" hose barb



- #69-21304 Ohmeda type vacuum adapter x D.I.S.S. vacuum female
- #69-20146 Ohmeda type oxygen adapter x D.I.S.S. oxygen female
- #69-20393 Ohmeda type medical air adapter x D.I.S.S. female
- #69-21886 Ohmeda type nitrous oxide adapter x D.I.S.S. nitrous oxide female
- #69-21305 Ohmeda type WAGD adapter x D.I.S.S. WAGD female
- #69-21306 Ohmeda carbon dioxide adapter x D.I.S.S. carbon dioxide female



- #69-20390 Ohmeda type oxygen adapter x 1/4 - 18 male
- #69-20142 Ohmeda type oxygen adapter x 1/8 - 27 male
- #69-20391 Ohmeda type medical air adapter x 1/4 - 18 male
- #69-20144 Ohmeda type medical air adapter x 1/8 - 27 male
- #69-20392 Ohmeda type vacuum adapter x 1/4 - 18 male
- #69-20145 Ohmeda type vacuum adapter x 1/8 - 27 male
- #69-20500 Ohmeda type WAGD adapter x 1/8 - 27 male
- #69-20143 Ohmeda type nitrous oxide adapter x 1/8 - 27 male
- #69-20761 Ohmeda carbon dioxide adapter x 1/8 - 27 male
- #69-21300 Ohmeda type WAGD adapter x 1/4 - 18 male
- #69-21301 Ohmeda carbon dioxide adapter x 1/4 - 18 male



- #69-21302 Ohmeda type WAGD adapter x D.I.S.S. WAGD male
- #69-30072 Ohmeda type oxygen adapter x D.I.S.S. oxygen male
- #69-30061 Ohmeda type vacuum adapter x D.I.S.S. vacuum male
- #69-20394 Ohmeda type nitrous oxide adapter x D.I.S.S. nitrous oxide male
- #69-30060 Ohmeda type medical air adapter x D.I.S.S. medical air male
- #69-21303 Ohmeda carbon dioxide adapter x D.I.S.S. carbon dioxide male



- #69-21528 Ohmeda type carbon dioxide check x 1/4 - 18 male
- #69-20147 Ohmeda type oxygen check x 1/4 - 18 male
- #69-20148 Ohmeda type nitrous oxide check x 1/4 - 18 male
- #69-20149 Ohmeda type air check x 1/4 - 18 male
- #69-20150 Ohmeda type vacuum check x - 1/4 - 18 male
- #69-20395 Ohmeda WAGD check x 1/4 - 18 male



- #69-20989 Ohmeda type oxygen coupler x 1/4" hose barb
- #69-20990 Ohmeda type nitrous oxide coupler x 1/4" hose barb
- #69-20991 Ohmeda type medical air coupler x 1/4" hose barb
- #69-20992 Ohmeda type vacuum coupler x 1/4" hose barb
- #69-21514 Ohmeda type vacuum coupler x 5/16" hose barb
- #69-20762 Ohmeda carbon dioxide coupler x 1/4" hose barb
- #69-21496 Ohmeda type WAGD coupler x 5/16" hose barb
- #69-21499 Ohmeda type WAGD coupler x 1/4" hose barb

Adapters/Fittings/ohmeda Style



- #69-30054 Ohmeda type oxygen male X (2) OHMEDA oxygen female checks
- #69-30055 Ohmeda type vacuum male x (2) OHMEDA vacuum female checks
- #69-30056 Ohmeda type nitrous oxide male x (2) OHMEDA nitrous oxide checks
- #69-30057 Ohmeda type medical air male (2) OHMEDA medical air checks



- #69-20270 Ohmeda type oxygen adapter x 1/8 female
- #69-20271 Ohmeda type nitrous oxide adapter x 1/8 female
- #69-20272 Ohmeda type medical air adapter x 1/8 female
- #69-20273 Ohmeda type vacuum adapter x 1/8 female
- #69-21297 Ohmeda type WAGD adapter x 1/8-27 female
- #69-21298 Ohmeda carbon dioxide adapter x 1/8-27 female



- #69-21524 Ohmeda WAGD coupler x 1/8 - 27 male
- #69-21337 Ohmeda oxygen coupler x 1/8 - 27 male
- #69-21338 Ohmeda nitrous oxide coupler x 1/8 - 27 male
- #69-21339 Ohmeda medical air coupler x 1/8 - 27 male
- #69-21340 Ohmeda vacuum coupler x 1/8 - 27 male
- #69-21526 Ohmeda carbon dioxide coupler x 1/8 - 27 male



- #69-21525 Ohmeda WAGD coupler x 1/8 - 27 female
- #69-21342 Ohmeda oxygen coupler x 1/8 - 27 female
- #69-21370 Ohmeda nitrous oxide coupler x 1/8 - 27 female
- #69-21371 Ohmeda medical air coupler x 1/8 - 27 female
- #69-21372 Ohmeda vacuum coupler x 1/8 - 27 female
- #69-21527 Ohmeda carbon dioxide coupler x 1/8 - 27 female



- #69-21373 Ohmeda oxygen coupler x D.I.S.S. oxygen male
- #69-21374 Ohmeda nitrous oxide coupler x D.I.S.S. nitrous oxide male
- #69-21375 Ohmeda medical air coupler x D.I.S.S. medical air male
- #69-21376 Ohmeda vacuum coupler x D.I.S.S. vacuum male



- #69-21863 Ohmeda oxygen coupler x D.I.S.S. oxygen female
- #69-21377 Ohmeda nitrous oxide coupler x D.I.S.S. nitrous oxide female
- #69-21864 Ohmeda medical air coupler x D.I.S.S. medical air female
- #69-21378 Ohmeda vacuum coupler x D.I.S.S. vacuum female

Adapters/Fittings/Oxequip Med*Star Style



- #69-20560 Oxequip WAGD Med*Star adapter x 1/8 - 27 male
- #69-20561 Oxequip WAGD Med*Star adapter x 1/4 - 18 male
- #69-20558 Oxequip carbon dioxide Med*Star adapter x 1/8 - 27 male
- #69-20559 Oxequip carbon dioxide Med*Star adapter x 1/4 - 18 male
- #69-20291 Oxequip oxygen Med*Star adapter x 1/8 - 27 male
- #69-20295 Oxequip oxygen Med*Star adapter x 1/4 - 18 male
- #69-20292 Oxequip nitrous oxide Med*Star adapter x 1/8 - 27 male
- #69-20296 Oxequip nitrous oxide Med*Star adapter x 1/4 - 18 male
- #69-20293 Oxequip medical air Med*Star adapter x 1/8 - 27 male
- #69-20297 Oxequip medical air Med*Star adapter x 1/4 - 18 male
- #69-20294 Oxequip vacuum Med*Star adapter x 1/8 - 27 male
- #69-20298 Oxequip vacuum Med*Star adapter x 1/4 - 18 male



- #69-20562 Oxequip WAGD Med*Star adapter x WAGD D.I.S.S. male
- #69-20301 Oxequip oxygen Med*Star adapter x oxygen D.I.S.S. male
- #69-20302 Oxequip nitrous oxide Med*Star adapter x nitrous oxide D.I.S.S. male
- #69-20303 Oxequip medical air Med*Star adapter x medical air D.I.S.S. male
- #69-20304 Oxequip vacuum Med*Star adapter x vacuum D.I.S.S. male
- #69-20563 Oxequip carbon dioxide Med*Star adapter x carbon dioxide D.I.S.S. male



- #69-21452 Oxequip oxygen Med*Star adapter x 1/8-27 female
- #69-21453 Oxequip nitrous oxide Med*Star adapter x 1/8-27 female
- #69-21454 Oxequip medical air Med*Star adapter x 1/8-27 female
- #69-21455 Oxequip vacuum Med*Star adapter x 1/8-27 female
- #69-21456 Oxequip carbon dioxide Med*Star adapter x 1/8-27 female
- #69-21457 Oxequip WAGD Med*Star adapter x 1/8-27 female



- #69-21458 Oxequip oxygen Med*Star x D.I.S.S. oxygen female
- #69-21459 Oxequip nitrous oxide med*adapter D.I.S.S. nitrous oxide female
- #69-21460 Oxequip medical air Med*Star adapter x D.I.S.S. medical air female
- #69-21461 Oxequip vacuum Med*Star adapter x D.I.S.S. vacuum female
- #69-21462 Oxequip carbon dioxide Med*Star adapter x D.I.S.S. carbon dioxide female
- #69-21463 Oxequip WAGD Med*Star adapter x D.I.S.S. WAGD female



- #69-20306 Oxequip oxygen Med*Star adapter x 90° oxygen D.I.S.S. male
- #69-20307 Oxequip nitrous oxide Med*Star adapter x 90° nitrous oxide D.I.S.S. male
- #69-20308 Oxequip medical air Med*Star adapter x 90° medical air D.I.S.S. male
- #69-20309 Oxequip vacuum Med*Star x 90° vacuum D.I.S.S. male



- #69-21824 Oxequip oxygen Med*Star adapter x 1/4" hose barb
- #69-21827 Oxequip nitrous oxide Med*Star adapter x 1/4" hose barb
- #69-21826 Oxequip medical air Med*Star adapter x 1/4" hose barb
- #69-21825 Oxequip vacuum Med*Star adapter x 1/4" hose barb
- #69-21341 Oxequip carbon dioxide Med*Star adapter x 1/4" hose barb
- #69-21828 Oxequip WAGD Med*Star adapter x 1/4" hose barb



- #69-21893 Oxequip oxygen Med*Star wye style duplex male adapter inlet x two Med*Star couplers
- #69-21465 Oxequip medical air Med*Star wye style duplex male adapter inlet x two Med*Star couplers
- #69-21894 Oxequip vacuum Med*Star wye style duplex male inlet x two Med*Star couplers

Note: Med*Star checks accept the original Oxequip twist style adapters, referred to as OES or 07 adapter - but original Oxequip OES checks or outlets do NOT accept Med*Star adapters.

Adapters/Fittings/Oxequip O.E.S. (07) Style



- #69-20151 Oxequip oxygen adapter x D.I.S.S. oxygen male
- #69-20152 Oxequip nitrous oxide adapter x D.I.S.S. nitrous oxide male
- #69-20153 Oxequip medical air adapter x D.I.S.S. medical air male
- #69-20154 Oxequip vacuum adapter x D.I.S.S. vacuum male
- #69-20951 Oxequip carbon dioxide adapter x D.I.S.S. carbon dioxide male
- #69-20952 Oxequip nitrogen adapter x D.I.S.S. nitrogen male
- #69-21316 Oxequip WAGD adapter x D.I.S.S. WAGD male



- #69-20155 Oxequip oxygen adapter x 1/8 - 27 male
- #69-21317 Oxequip oxygen adapter x 1/4 - 18 male
- #69-20156 Oxequip nitrous oxide adapter x 1/8 - 27 male
- #69-21319 Oxequip nitrous oxide adapter x 1/4 - 18 male
- #69-20157 Oxequip medical air adapter x 1/8 - 27 male
- #69-21318 Oxequip medical air adapter x 1/4 - 18 male
- #69-20158 Oxequip vacuum adapter x 1/8 - 27 male
- #69-21320 Oxequip vacuum adapter x 1/4 - 18 male
- #69-20953 Oxequip carbon dioxide adapter x 1/8 - 27 male
- #69-20954 Oxequip carbon dioxide adapter x 1/4 - 18 male
- #69-20955 Oxequip nitrogen adapter x 1/8 - 27 male
- #69-20956 Oxequip nitrogen adapter x 1/4 - 18 male



- #69-20284 Oxequip oxygen adapter x D.I.S.S. oxygen female
- #69-21321 Oxequip nitrous oxide adapter x D.I.S.S. nitrous oxide female
- #69-21322 Oxequip medical air adapter x D.I.S.S. medical air female
- #69-21323 Oxequip vacuum adapter x D.I.S.S. vacuum female
- #69-21324 Oxequip carbon dioxide adapter x D.I.S.S. carbon dioxide female
- #69-21325 Oxequip nitrogen adapter x D.I.S.S. nitrogen female
- #69-21326 Oxequip WAGD adapter x D.I.S.S. WAGD female



- #69-20285 Oxequip oxygen adapter x 90° elbow D.I.S.S. oxygen male
- #69-20286 Oxequip nitrous oxide adapter x 90° elbow D.I.S.S. nitrous oxide male
- #69-20287 Oxequip medical air adapter x 90° elbow D.I.S.S. medical air male
- #69-20288 Oxequip vacuum adapter x 90° elbow D.I.S.S. vacuum male



- #69-20419 Oxequip oxygen check x D.I.S.S. oxygen male
- #69-20420 Oxequip nitrous oxide check x D.I.S.S. nitrous oxide male
- #69-20421 Oxequip medical air check x D.I.S.S. vacuum male
- #69-20422 Oxequip vacuum check x D.I.S.S. vacuum male
- #69-20159 Oxequip oxygen check x 1/4 -18 male
- #69-20161 Oxequip medical air check x 1/4 - 18 male
- #69-20162 Oxequip vacuum check x 1/4 - 18 male
- #69-20160 Oxequip nitrous oxide check x 1/4 - 18 male



- #69-21333 Oxequip oxygen adapter x 1/4" hose barb
- #69-21334 Oxequip nitrous oxide adapter x 1/4" hose barb
- #69-21335 Oxequip medical air adapter x 1/4" hose barb
- #69-21336 Oxequip vacuum adapter x 1/4" hose barb



- #69-21327 Oxequip oxygen adapter x 1/8 - 27 female
- #69-21328 Oxequip nitrous oxide adapter x 1/8 - 27 female
- #69-21329 Oxequip medical air adapter x 1/8 - 27 female
- #69-21330 Oxequip vacuum adapter x 1/8 - 27 female
- #69-21331 Oxequip carbon dioxide adapter x 1/8 - 27 female
- #69-21332 Oxequip nitrogen adapter x 1/8 - 27 female

Note: Med*Star checks accept the original Oxequip twist style adapters, referred to as OES or 07 adapter - but original Oxequip OES checks or outlets do NOT accept Med*Star adapters.

Adapters/Fittings/Puritan-Bennett



- #69-20163 Puritan oxygen adapter x 1/4" hose barb
- #69-20164 Puritan nitrous oxide adapter x 1/4" hose barb
- #69-20165 Puritan medical air adapter x 1/4" hose barb
- #69-20166 Puritan vacuum adapter x 1/4" hose barb
- #69-21520 Puritan vacuum adapter x 5/16" hose barb
- #69-20396 Puritan WAGD adapter x 5/16" hose barb
- #69-21481 Puritan WAGD adapter x 1/4" hose barb



- #69-20167 Puritan type oxygen check x 1/4" hose barb
- #69-20168 Puritan type nitrous oxide check x 1/4" hose barb
- #69-20169 Puritan type medical air check x 1/4" hose barb
- #69-20170 Puritan type vacuum check x 1/4" hose barb
- #69-21518 Puritan type vacuum check x 5/16" hose barb
- #69-20397 Puritan type WAGD check x 1/4" hose barb
- #69-21519 Puritan type WAGD check x 5/16" hose barb



- #69-20355 Puritan oxygen adapter x 1/4 - 18 male
- #69-20171 Puritan oxygen adapter x 1/8 - 27 male
- #69-21307 Puritan nitrous oxide adapter x 1/8 - 27 male
- #69-20172 Puritan nitrous oxide adapter x 1/4 - 18 male
- #69-20173 Puritan medical air adapter x 1/4 - 18 male
- #69-21308 Puritan medical air adapter x 1/8 - 27 male
- #69-20174 Puritan vacuum adapter x 1/8 - 27 male
- #69-20356 Puritan vacuum adapter x 1/4 - 18 male
- #69-21248 Puritan WAGD adapter x 1/4 - 18 male
- #69-21545 Puritan WAGD adapter x 1/8 - 27 male



- #69-21309 Puritan nitrous oxide adapter x 1/8 - 27 female
- #69-21310 Puritan medical air adapter x 1/8 - 27 female
- #69-20274 Puritan oxygen adapter x 1/8 - 27 female
- #69-21311 Puritan vacuum adapter x 1/8 - 27 female
- #69-21546 Puritan WAGD adapter x 1/8 - 27 female



- #69-20275 Puritan oxygen adapter x oxygen D.I.S.S. male
- #69-20276 Puritan vacuum adapter x vacuum D.I.S.S. male
- #69-20398 Puritan medical air adapter x medical air D.I.S.S. male
- #69-20621 Puritan nitrous oxide adapter x nitrous oxide D.I.S.S. male
- #69-21544 Puritan WAGD adapter x D.I.S.S. male



- #69-20277 Puritan type oxygen check x 1/4 - 18 male
- #69-20278 Puritan type nitrous oxide check x 1/4 - 18 male
- #69-20279 Puritan type medical air check x 1/4 - 18 male
- #69-20280 Puritan type vacuum check 1/4 - 18 male
- #69-21250 Puritan type WAGD check x 1/4 - 18 male



- #69-20281 Puritan oxygen duplex adapter, Puritan inlet stem and (2) Puritan check outlets
- #69-20282 Puritan medical air duplex adapter, Puritan inlet stem and (2) Puritan check outlets
- #69-20357 Puritan vacuum duplex adapter, Puritan inlet stem and (2) Puritan check outlets



- #69-21312 Puritan oxygen adapter x D.I.S.S. oxygen female
- #69-21313 Puritan nitrous oxide adapter x D.I.S.S. nitrous oxide female
- #69-21314 Puritan medical air adapter x D.I.S.S. medical air female
- #69-21315 Puritan vacuum adapter x D.I.S.S. vacuum female
- #69-21547 Puritan WAGD adapter x D.I.S.S. WAGD female

Adapters/Fittings/Schrader Style

NON-Swivel DESIGN - Used where positioning is critical - Example: Flowmeters and Suction Regulators



- #69-20133 Schrader oxygen adapter x 1/8-27 male
- #69-20410 Schrader nitrous oxide adapter x 1/8 - 27 male
- #69-20414 Schrader medical air adapter x 1/8 - 27 male*
- #69-20412 Schrader vacuum adapter x 1/8 - 27 male



- #69-21387 Schrader oxygen adapter x 1/8-27 female
- #69-21388 Schrader nitrous oxide adapter x 1/8-27 female
- #69-21389 Schrader medical air adapter x 1/8-27 female*
- #69-21390 Schrader vacuum adapter x 1/8-27 female



- #69-20358 Schrader vacuum adapter x 1/4 - 18 male
- #69-20406 Schrader oxygen adapter x 1/4 - 18 male
- #69-20359 Schrader medical air adapter x 1/4 - 18 male*
- #69-20411 Schrader nitrous oxide adapter x 1/4 - 18 male



- #69-21413 Schrader oxygen adapter x 1/4" hose barb
- #69-21414 Schrader nitrous oxide adapter x 1/4" hose barb
- #69-21415 Schrader medical air adapter x 1/4" male hose barb*
- #69-21416 Schrader vacuum adapter x 1/4" hose barb
- #69-21515 Schrader vacuum adapter x 5/16" hose barb



- #69-20132 Schrader oxygen adapter x D.I.S.S. oxygen male
- #69-21236 Schrader nitrous oxide adapter x D.I.S.S. nitrous oxide male
- #69-21237 Schrader medical air adapter x D.I.S.S. medical air male
- #69-21238 Schrader vacuum adapter x D.I.S.S. vacuum male



- #69-20416 Schrader oxygen adapter x D.I.S.S. oxygen female
- #69-21417 Schrader nitrous oxide adapter x D.I.S.S. nitrous oxide female
- #69-21418 Schrader medical air adapter x D.I.S.S. medical air female
- #69-21419 Schrader vacuum adapter x D.I.S.S. vacuum female



- #83-900-0212 Nitrogen Hazard Warning Tag

Important Industry Update:

For years Schrader adapters and check units were not manufactured with nitrogen-specific indexing. In 2004, Schrader released such fittings which are now recommended for all nitrogen-driven products in place of the medical air version adopted by the industry for decades. Please use care when ordering.

*If used for nitrogen, specific potential hazard warnings should be affixed to device indicating such.

Adapters/Fittings/Schrader Style

SWIVEL DESIGN - USED ON HOSES AND OTHER ASSEMBLIES WHERE POSITIONING IS NOT REQUIRED



- #69-20408 Schrader oxygen adapter x 1/8-27 male
- #69-20409 Schrader nitrous oxide adapter x 1/8-27 male
- #69-20413 Schrader medical air adapter x 1/8-27 male*
- #69-20404 Schrader vacuum adapter x 1/8-27 male




- #69-21420 Schrader oxygen adapter x 1/8-27 female
- #69-21421 Schrader nitrous oxide adapter x 1/8-27 female
- #69-21422 Schrader medical air adapter x 1/8-27 female*
- #69-21423 Schrader vacuum adapter x 1/8-27 female



- #69-21535 Schrader nitrogen adapter x 1/4-18 male
- #69-21255 Schrader oxygen adapter x 1/4-18 male
- #69-21495 Schrader nitrous oxide adapter x 1/4-18 male
- #69-20415 Schrader medical air adapter x 1/4-18 male*
- #69-20405 Schrader vacuum adapter x 1/4-18 male



- #69-20128 Schrader oxygen adapter x 1/4" hose barb
- #69-20129 Schrader nitrous oxide adapter x 1/4" hose barb
- #69-20130 Schrader medical air adapter x 1/4" hose barb*
- #69-20131 Schrader vacuum adapter x 1/4" hose barb
- #69-21516 Schrader vacuum adapter x 5/16" hose barb
- #69-21534 Schrader nitrogen adapter x 1/4" hose barb



- #69-21424 Schrader oxygen adapter x D.I.S.S. oxygen male
- #69-21425 Schrader nitrous oxide adapter x D.I.S.S. nitrous oxide male
- #69-21426 Schrader medical air adapter x D.I.S.S. medical air male
- #69-21427 Schrader vacuum adapter x D.I.S.S. vacuum male



- #69-21428 Schrader oxygen adapter x D.I.S.S. oxygen female
- #69-21429 Schrader nitrous oxide adapter x D.I.S.S. nitrous oxide female
- #69-21430 Schrader medical air adapter x D.I.S.S. medical air female
- #69-21431 Schrader vacuum adapter x D.I.S.S. vacuum female



- #83-900-0212 Nitrogen Hazard Warning Tag

Important Industry Update:

For years Schrader adapters and check units were not manufactured with nitrogen-specific indexing. In 2004, Schrader released such fittings which are now recommended for all nitrogen-driven products in place of the medical air version adopted by the industry for decades. Please use care when ordering.

*If used for nitrogen, specific potential hazard warnings should be affixed to device indicating such.

Adapters/Fittings/Schrader Style

	#69-20134	Schrader oxygen check 1/4-18 male
	#69-20135	Schrader nitrous oxide check x 1/4-18 male
	#69-20136	Schrader medical air check x 1/4-18 male (also used for nitrogen)*
	#69-20137	Schrader vacuum check x 1/4-18 male
	#69-21533	Schrader nitrogen check x 1/4-18 male
	#69-21432	Schrader oxygen check x 1/8-27 male
	#69-21433	Schrader nitrous oxide check x 1/8-27 male*
	#69-21434	Schrader medical air check x 1/8-27 male
	#69-21435	Schrader vacuum check x 1/8-27 male
	#69-21436	Schrader oxygen check x 1/4" hose barb
	#69-21437	Schrader nitrous oxide check x 1/4" hose barb
	#69-21438	Schrader medical air check x 1/4" hose barb*
	#69-21439	Schrader vacuum check x 1/4" hose barb
	#69-21517	Schrader vacuum check x 5/16" hose barb
	#69-21549	Schrader nitrogen check x 1/4" hose barb
	#69-21440	Schrader oxygen check x 1/8-27 female
	#69-21441	Schrader nitrous oxide check x 1/8-27 female
	#69-21442	Schrader medical air check x 1/8-27 female*
	#69-21443	Schrader vacuum check x 1/8-27 female
	#69-21444	Schrader oxygen check x D.I.S.S. oxygen male
	#69-21445	Schrader nitrous oxide check x D.I.S.S. nitrous oxide male
	#69-21446	Schrader medical air check x D.I.S.S. medical air male
	#69-21447	Schrader vacuum check x D.I.S.S. vacuum male
	#69-21448	Schrader oxygen check x D.I.S.S. oxygen female
	#69-21449	Schrader nitrous oxide check x D.I.S.S. nitrous oxide female
	#69-21450	Schrader medical air check x D.I.S.S. medical air female
	#69-21451	Schrader vacuum check x D.I.S.S. vacuum female
	#69-21536	Schrader nitrogen check x D.I.S.S. nitrogen female
	#69-20283	Schrader oxygen duplex adapter, Schrader inlet stem and two Schrader check outlets
	#69-20362	Schrader medical air duplex adapter, Schrader inlet stem and two Schrader check outlets*
	#69-20418	Schrader vacuum duplex adapter, Schrader inlet and two Schrader check outlets
	#69-21548	Schrader nitrogen duplex adapter, Schrader inlet stem and two Schrader nitrogen check outlets

*If used for nitrogen, specific hazard warnings should be affixed to device indicating such, refer to Hazard Warning label #83-900-0212.

Adapters/Fittings/D.I.S.S.



- #69-20185 Oxygen hex nut
- #69-20186 Nitrous oxide hex nut
- #69-20187 Medical air hex nut
- #69-20188 Vacuum hex nut
- #69-20324 Nitrogen hex nut
- #69-20434 Helium hex nut
- #69-20435 Carbon dioxide hex nut
- #69-20959 WAGD hex nut
- #69-21261 Heliox hex nut, helium > 80.5% (CGA #1180A)
- #69-21556 Instrument air hex nut



- #69-30200 Oxygen hand nut
- #69-30201 Nitrous oxide hand nut
- #69-30202 Medical air hand nut
- #69-21581 Instrument air hand nut
- #69-30203 Vacuum hand nut
- #69-30204 Nitrogen/ Helium hand nut
- #69-20437 Carbon dioxide hand nut
- #69-20566 WAGD hand nut




- #69-30300 Retaining ring for nitrous oxide, nitrogen, vacuum, medical air carbon dioxide and helium nipple
- #69-30302 Retaining ring for WAGD nipple



- #69-20218 Nipple, oxygen x 1/8 - 27 male
- #69-20219 Nipple, nitrous oxide x 1/8 - 27 male
- #69-20220 Nipple, nitrous oxide x 1/4 - 18 male
- #69-20221 Nipple, medical air x 1/8 - 27 male
- #69-20222 Nipple, medical air x 1/4 - 18 male
- #69-20223 Nipple, vacuum x 1/8 - 27 male
- #69-20224 Nipple, vacuum x 1/4 - 18 male
- #69-20322 Nipple, nitrogen x 1/8 - 27 male
- #69-20323 Nipple, nitrogen x 1/4 - 18 male
- #69-20995 Nipple, carbon dioxide x 1/8 - 27 male
- #69-20996 Nipple, carbon dioxide x 1/4 - 18 male
- #69-21491 Nipple, WAGD x 1/8 - 27 male
- #69-21553 Nipple, instrument air x 1/8 - 27 male
- #69-21554 Nipple, instrument air x 1/4 - 18 male



- #69-20432 Oxygen hose barb for 3/16" I.D. hose
- #69-20181 Oxygen hose barb for 1/4" I.D. hose
- #69-20182 Nitrous oxide hose barb for 1/4" I.D. hose
- #69-20183 Medical air hose barb for 1/4" I.D. hose
- #69-20184 Vacuum hose barb for 1/4" I.D. hose
- #69-21501 Vacuum hose barb for 5/16" I.D. hose
- #69-30301 Nitrogen hose barb for 1/4" I.D. hose
- #69-20430 Helium hose barb for 1/4" I.D. hose
- #69-20431 Carbon dioxide hose barb for 1/4" I.D. hose
- #69-20960 WAGD hose barb for 1/4" I.D. hose
- #69-21500 WAGD hose barb for 5/16" I.D. hose
- #69-21262 Heliox hose barb for 1/4" I.D. hose, helium > 80.5% (CGA #1180A)
- #69-21555 Instrument air hose barb for 1/4" I.D. hose



IT PIVOTS!

- #69-21577 Medical Air D.I.S.S. female x 90° pivoting elbow with 1/4" hose barb
- #69-21578 Oxygen D.I.S.S. female x 90° pivoting elbow with 1/4" hose barb
- #69-21579 Nitrous Oxide D.I.S.S. female x 90° pivoting elbow with 1/4" hose barb
- #69-21580 Vacuum D.I.S.S. female x 90° pivoting elbow with 1/4" hose barb

Adapters/Fittings/D.I.S.S./Female




Chain with hook



#69-21240	Oxygen handnut and nipple combination for 1/4" I.D. hose
#69-20439	Nitrous oxide handnut and nipple combination for 1/4" I.D. hose
#69-20440	Medical air handnut and nipple combination for 1/4" I.D. hose
#69-21582	Instrument air handnut and nipple for 1/4" I.D. hose
#69-20441	Vacuum handnut and nipple combination for 1/4" I.D. hose
#69-21506	Vacuum handnut and nipple combination for 5/16" I.D. hose
#69-20442	Nitrogen handnut and nipple combination for 1/4" I.D. hose
#69-20701	Helium handnut and nipple combination for 1/4" I.D. hose
#69-20443	Carbon dioxide handnut nipple combination for 1/4" I.D. hose
#69-21245	WAGD handnut and nipple combination for 1/4" I.D. hose
#69-21507	WAGD handnut and nipple combination for 5/16" I.D. hose
#69-21852	95% O ₂ , 5% CO ₂ handnut and nipple combination for 1/4" I.D. hose
#69-21343	Oxygen nut and nipple combination x 1/8 - 27 male
#69-21344	Nitrous oxide and nipple combination x 1/8 - 27 male
#69-21345	Medical air nut and nipple combination x 1/8 - 27 male
#69-21346	Vacuum nut and nipple combination x 1/8 - 27 male
#69-21347	Carbon dioxide nut and nipple combination x 1/8 - 27 male
#69-21348	WAGD nut and nipple combination x 1/8 - 27 male
#69-21349	Nitrogen nut and nipple combination x 1/8 - 27 male
#69-21350	95% O ₂ , 5% CO ₂ nut and nipple combination x 1/8 - 27 male
#69-21550	Instrument air nut and nipple combination x 1/8 - 27 male
#69-21551	Instrument air nut and nipple combination x 1/4 - 18 male
#69-20325	Oxygen nut nipple combination for 1/4" I.D. hose
#69-20326	Nitrous oxide nut and nipple combination for 1/4" I.D. hose
#69-20327	Medical air nut and nipple combination for 1/4" I.D. hose
#69-20328	Vacuum nut and nipple combination for 1/4" I.D. hose
#69-21504	Vacuum nut and nipple combination for 5/16" I.D. hose
#69-20329	Nitrogen nut and nipple combination for 1/4" I.D. hose
#69-20700	Helium nut and nipple combination for 1/4" I.D. hose
#69-20438	Carbon dioxide nut and nipple combination for 1/4" I.D. hose
#69-21505	WAGD nut and nipple combination for 5/16" I.D. hose
#69-21246	WAGD nut and nipple combination for 1/4" I.D. hose
#69-21815	95% O ₂ , 5% CO ₂ nut and nipple combination for 1/4" I.D. hose
#69-21816	Special gas mixture nut and nipple combination for 1/4" I.D. hose (CGA 1020A)
#69-21552	Instrument air nut and nipple combination for 1/4" I.D. hose
#69-21264	Oxygen nut and nipple combination x tapered hose barb, 3/16" to 3/8", with 5" chain and hook
#69-20535	Oxygen D.I.S.S. plug and chain assembly
#69-20993	WAGD D.I.S.S. plug and chain assembly
#69-21247	Common D.I.S.S. plug and chain assembly for nitrous oxide, medical air, vacuum, instrument air and nitrogen
#69-21359	Nitrous oxide handnut and nipple combination x 1/8 - 27 female
#69-21360	Medical air handnut and nipple combination x 1/8 - 27 female
#69-21362	Carbon dioxide handnut and nipple combination x 1/8 - 27 female
#69-21363	WAGD handnut and nipple combination x 1/8 - 27 female
#69-21364	Nitrogen handnut and nipple combination x 1/8 - 27 female
#69-21365	95% O ₂ , 5% CO ₂ handnut and nipple combination x 1/8 - 27 female
#69-20727	Oxygen D.I.S.S. plastic dust cap w/chain and hook
#69-20728	Oxygen D.I.S.S. plastic dust cap w/chain and eyelet
#69-20729	All other gases D.I.S.S. plastic dust cap w/chain and hook
#69-20730	All other gases D.I.S.S. plastic dust cap w/chain and eyelet
#69-21351	Oxygen handnut and nipple combination x 1/8 - 27 male
#69-21352	Nitrous oxide handnut and nipple combination x 1/8 - 27 male
#69-21353	Medical air handnut and nipple combination x 1/8 - 27 male
#69-21354	Vacuum handnut and nipple combination x 1/8 - 27 male
#69-21355	Carbon dioxide handnut and nipple combination x 1/8 - 27 male
#69-21356	WAGD handnut and nipple combination x 1/8 - 27 male
#69-21357	Nitrogen handnut and nipple combination x 1/8 - 27 male
#69-21358	95% O ₂ , 5% CO ₂ handnut and nipple combination x 1/8 - 27 male

Adapters/Fittings/D.I.S.S./Male



- #69-21575 Instrument Air D.I.S.S. check valve w/red hand grip x 1/4" hose barb
- #69-21509 WAGD D.I.S.S. check valve w/purple hand grip x 5/16" hose barb
- #69-21490 WAGD D.I.S.S. check valve w/purple hand grip x 1/4" hose barb
- #69-20879 Oxygen D.I.S.S. check valve w/green hand grip x 1/4" hose barb
- #69-20880 Nitrous oxide D.I.S.S. check valve w/blue hand grip x 1/4" hose barb
- #69-20881 Medical air D.I.S.S. check valve w/yellow hand grip x 1/4" hose barb
- #69-20882 Vacuum D.I.S.S. check valve w/white hand grip x 1/4" hose barb
- #69-21508 Vacuum D.I.S.S. check valve w/white hand grip x 5/16" hose barb
- #69-20883 Nitrogen D.I.S.S. check valve w/black hand grip x 1/4" hose barb
- #69-20998 Carbon dioxide D.I.S.S. check valve w/grey hand grip x 1/4" hose barb




- #69-20884 Oxygen D.I.S.S. male adapter w/green hand grip x 1/4" hose barb
- #69-20885 Nitrous oxide D.I.S.S. male adapter w/blue hand grip x 1/4" hose barb
- #69-20886 Medical air D.I.S.S. male adapter w/yellow hand grip x 1/4" hose barb
- #69-20888 Nitrogen D.I.S.S. male adapter w/black hand grip x 1/4" hose barb
- #69-20994 Carbon dioxide D.I.S.S. male adapter w/grey hand grip x 1/4" hose barb
- #69-20887 Vacuum D.I.S.S. male adapter w/white hand grip x 1/4" hose barb
- #69-21503 Vacuum D.I.S.S. male adapter w/white hand grip x 5/16" hose barb
- #69-21466 WAGD D.I.S.S. male adapter w/purple hand grip x 1/4" hose barb
- #69-21502 WAGD D.I.S.S. male adapter w/purple hand grip x 5/16" hose barb




- #69-20201 Oxygen adapter, D.I.S.S. male x 1/4" hose barb
- #69-30101 Nitrous oxide adapter, D.I.S.S. male x 1/4" hose barb
- #69-30102 Medical air adapter, D.I.S.S. male x 1/4" hose barb
- #69-30103 Vacuum adapter, D.I.S.S. male x 1/4" hose barb
- #69-20330 Nitrogen adapter, D.I.S.S. male x 1/4" hose barb



- #69-20175 Adapter, oxygen D.I.S.S. male x oxygen D.I.S.S. male
- #69-20204 Adapter, oxygen D.I.S.S. male x 1/8 - 27 male
- #69-20177 Adapter, oxygen D.I.S.S. male x 1/4 - 18 male



- #69-20230 Adapter, nitrous oxide D.I.S.S. male x 1/8 - 27 male
- #69-20231 Adapter, nitrous oxide D.I.S.S. male x 1/4 - 18 male
- #69-20232 Adapter, medical air D.I.S.S. male x 1/8 - 27 male
- #69-20233 Adapter, medical air D.I.S.S. male x 1/4 - 18 male
- #69-20234 Adapter, vacuum D.I.S.S. male x 1/8 - 27 male
- #69-20235 Adapter, vacuum D.I.S.S. male x 1/4 - 18 male
- #69-20319 Adapter, nitrogen D.I.S.S. male x 1/8 - 27 male
- #69-20320 Adapter, nitrogen D.I.S.S. male x 1/4 - 18 male
- #69-20428 Adapter, helium D.I.S.S. male x 1/4 - 18 male
- #69-20429 Adapter, carbon dioxide D.I.S.S. male x 1/4 - 18 male
- #69-21571 Adapter carbon dioxide D.I.S.S. male x 1/8 - 27 male
- #69-20961 Adapter, WAGD D.I.S.S. male x 1/4 - 18 male
- #69-21557 Adapter, instrument air D.I.S.S. male x 1/8 - 27 male
- #69-21558 Adapter, instrument air D.I.S.S. male x 1/4 - 18 male



- #69-20210 Adapter, oxygen D.I.S.S. male x 1/8 - 27 female
- #69-20211 Adapter, oxygen D.I.S.S. male x 1/4 - 18 female



- #69-20236 Adapter, nitrous oxide D.I.S.S. male x 1/8 - 27 female
- #69-20237 Adapter, medical air D.I.S.S. male x 1/8 - 27 female
- #69-20321 Adapter, nitrogen D.I.S.S. male x 1/8 - 27 female
- #69-20238 Adapter, vacuum D.I.S.S. male x 1/8 - 27 female



- #69-21216 Coupler grip 1/8 - 27 NPT female x 1/8 - 27 NPT male - Oxygen
- #69-21217 Coupler grip 1/8 - 27 NPT female x 1/8 - 27 NPT male - Medical Air
- #69-21218 Coupler grip 1/8 - 27 NPT female x 1/8 - 27 NPT male - Vacuum
- #69-21219 Coupler grip 1/8 - 27 NPT female x 1/8 - 27 NPT Male - Nitrous Oxide

Adapters/Fittings/D.I.S.S./Female X D.I.S.S. Male



- #69-20750 - Swivel adapter, 90°, D.I.S.S. oxygen female X 1/8-27 female (not shown)
- #69-20751 - Swivel adapter, 90° D.I.S.S. oxygen female X 1/4-18 female (not shown)
- #69-20638 - Swivel adapter, 90°, D.I.S.S. oxygen female X 1/8-27 male - 1 1/2" long
- #69-20639 - Swivel adapter, 90°, D.I.S.S. oxygen female X 1/4-18 male - 1 1/2" long



- #69-20640 - Swivel adapter, D.I.S.S. oxygen female x 1/8-27 female
- #69-20641 - Swivel adapter, D.I.S.S. oxygen female x 1/4-18 female



- #69-20642 - Swivel adapter, D.I.S.S. oxygen female x 1/8-27 male
- #69-20643 - Swivel adapter, D.I.S.S. oxygen female x 1/4-18 male




- #69-20178 - Coupling - D.I.S.S. female oxygen x D.I.S.S. female oxygen



- #69-21805 - Swivel elbow, D.I.S.S. female oxygen x D.I.S.S. female oxygen



- #69-20636 - Swivel elbow, D.I.S.S. female oxygen x D.I.S.S. male oxygen (3/4" long)
- #69-20217 - Swivel elbow, D.I.S.S. female oxygen x D.I.S.S. male oxygen (1 1/2" long)
- #69-20637 - Swivel elbow, D.I.S.S. female oxygen x D.I.S.S. male oxygen (2" long)
- #69-20353 - Swivel elbow, D.I.S.S. female medical air x D.I.S.S. male medical air
- #69-21879 - Swivel elbow, D.I.S.S. female vacuum x D.I.S.S. male vacuum



- #69-21260 - Adapter, tee 1/8-27 NPT female inlet x two D.I.S.S. female oxygen outlets



- #69-20644 - Oxygen ... D.I.S.S. 1240 Triple Outlet Rail Block w/check valve outlets
- #69-20645 - Nitrous Oxide ... D.I.S.S. 1040-A Triple Outlet Rail Block w/ D.I.S.S. male outlets
- #69-20646 - Air ... D.I.S.S. 1160-A Triple Outlet Rail Block w/ D.I.S.S. male outlets
- #69-20647 - Nitrogen ... D.I.S.S. 1120-A Triple Outlet Rail Block w/ D.I.S.S. male outlets
- #69-20648 - Rail block less fittings ... 1/8" NPT female inlet ... (3) 1/4" NPT female outlets w/ mounting bracket
- #69-20649 - Female bracket for wall mounting

Adapters/Fittings/D.I.S.S.



- #69-20193 Oxygen duplex adapter w/o check valves
- #69-20194 Oxygen duplex adapter w/check valves
- #69-20195 Medical air duplex adapter w/o check valves
- #69-20196 Medical air duplex adapter w/check valves
- #69-20197 Vacuum duplex adapter w/o check valves
- #69-20198 Vacuum duplex adapter w/check valves
- #69-20619 Nitrogen duplex adapter w/o check valves
- #69-21806 Nitrogen duplex adapter w/check valves
- #69-20999 Nitrous oxide duplex adapter w/o check valves
- #69-21200 Nitrous oxide duplex adapter w/check valves
- #69-21537 Carbon dioxide adapter w/check valves
- #69-21559 Carbon dioxide duplex adapter w/check valves



- #69-20586 Replacement O-Ring for D.I.S.S. oxygen nipple 10/pk
- #69-20587 Replacement O-Ring for D.I.S.S. nitrous oxide nipple 10/pk
- #69-20588 Replacement O-Ring for D.I.S.S. nitrogen nipple 10/pk
- #69-20589 Replacement O-Ring for D.I.S.S. vacuum nipple 10/pk
- #69-20590 Replacement O-Ring for D.I.S.S. medical air nipple 10/pk
- #69-21489 Replacement O-Ring for D.I.S.S. WAGD nipple 10/pk



- #69-20433 Oxygen hose barb, tapered for 3/16" to 3/8" I.D. hose



- #69-20565 Oxygen wing nut



- #69-20921 Vacuum D.I.S.S. extension assembly. This D.I.S.S. female to D.I.S.S. male assembly is used between D.I.S.S. outlet of suction regulator and D.I.S.S. inlet type disposable collection unit to unit to allow replacement of disposable collection unit without the need to remove the suction regulator from the wall. Also serves to overcome space problems. Fits "fork" bracket adapters.



- #69-20922 Vacuum D.I.S.S. to stem assembly. This D.I.S.S. female to serrated stem assembly allows the user to mount a suction regulator on the "eye" hook of a floor mount holder. This permits use of the suction regulator on the opening room floor when stem inlet type disposable collection units are preferred. Fits "fork" bracket adapters.

Adapters/Fittings/Check Valves





DEMAND CHECK VALVES

50 PSI Working Pressure

Demand Check Valve Units are designed for use in D.I.S.S. outlet systems where immediate gas supply (or suction) is desired on demand. Positive shut off is achieved with disconnection of inlet nut and nipple. Demand Check Valve Units are to serve secondary equipment downstream from standard service outlets.

Full free extension of poppet is required for unobstructed full flow pressure when demand check is connected to D.I.S.S. hose nut and nipple. If proper clearance is not possible, then Extension Adapters are necessary. After installation, test for free full flow of pressure prior to operation.

#69-20628 - Typical Operation



#69-21256 - Typical Operation



- #69-21576 Check valve, instrument air D.I.S.S. male x 1/4 - 18 male
- #69-21256 Check valve, WAGD D.I.S.S. male x 1/4 - 18 male
- #69-20239 Check valve, oxygen D.I.S.S. male x 1/8 - 27 male
- #69-20240 Check valve, oxygen D.I.S.S. male x 1/4 - 18 male
- #69-20241 Check valve, nitrous oxide D.I.S.S. male x 1/8 - 27 male
- #69-20242 Check valve, nitrous oxide D.I.S.S. male 1/4 - 18 male
- #69-20243 Check valve, medical air D.I.S.S. male x 1/8 - 27 male
- #69-20244 Check valve, medical air D.I.S.S. male x 1/4 - 18 male
- #69-20245 Check valve, vacuum D.I.S.S. male x 1/4 - 18 male
- #69-20317 Check valve, nitrogen D.I.S.S. male x 1/8 - 27 male
- #69-20318 Check valve, nitrogen D.I.S.S. male x 1/4 - 18 male
- #69-20997 Check valve, carbon dioxide D.I.S.S. male 1/4 - 18 male
- #69-21259 Check valve, carbon dioxide D.I.S.S. male x 1/8 - 27 male
- #69-21583 Check valve, helium-oxygen (Heliox) He>80% D.I.S.S. male x 1/8 - 27 male
- #69-21584 Check valve, helium-oxygen (Heliox) He>80% D.I.S.S. male x 1/4 - 18 male



- #69-20246 Check extension adapter, 1/8 - 27 male x 1/8 - 27 female
- #69-20247 Check extension adapter, 1/4 - 18 male x 1/8 - 27 female
- #69-20248 Check extension adapter, 1/4 - 18 male x 1/4 - 18 female




- #69-20628 Check valve - one way, oxygen D.I.S.S. male x 1/8 - 27 male
- #69-20629 Check valve - one way, oxygen D.I.S.S. male x 1/4 - 18 male
- #69-20630 Check valve - one way, nitrous oxide D.I.S.S. male x 1/8 - 27 male
- #69-20631 Check valve - one way, nitrous oxide D.I.S.S. male x 1/4 - 18 male
- #69-20632 Check valve - one way, medical air D.I.S.S. male x 1/8 - 27 male
- #69-20633 Check valve - one way, medical air D.I.S.S. male x 1/4 - 18 male
- #69-20635 Check valve - one way, nitrogen D.I.S.S. male x 1/4 - 18 male

CHECK VALVES - CHECK REVERSE FLOW


Design: These valves are designed to allow free flow of gas through supply lines and to check any reverse flow. Check valves are opened for full flow at only 4 oz. of gas pressure. If reverse flow occurs - the valve closes instantly with total shut-off at only 10 oz. of back pressure. Check valves are nickel plated.

Performance: Check valves have minimum flow volume of 1000 cu. ft/hr. at average of 50 P.S.I. Recommended working pressure not to exceed 125 P.S.I.


Safety: Valves should be inspected regularly and replaced when tests indicate faulty operation. Careless usage, dirt or abuse can shorten service life. All check valves are inspected before packaging.



#69-20553




#69-20554




#69-20555



#69-20556



#69-20557



#69-20622

CGA CONNECTION NO. 1240 OXYGEN

- #69-20553 Oxygen 9/16" - 18 R.H Stationary Female Nut to Male 9/16" - 18 Checks from Female to Male
- #69-20554 Oxygen 9/16" - 18 R.H Male to 9/16" - 18 Female Stationary Nut - Checks from Male to Female
- #69-20555 Oxygen 9/16" - 18 R.H Male to 1/4" NPT Male - Checks from D.I.S.S. side to 1/4" NPT Male
- #69-20556 All Gases 1/4" NPT Male to 1/4" NPT Male
- #69-20557 All Gases 1/4" NPT Female to 1/4" NPT Male - Checks from Female to Male
- #69-20622 All Gases 1/4" NPT Female to 1/4" NPT Male - Checks from Male to Female
- #69-21538 All Gases 1/4" NPT Female to 1/4" NPT Female (not shown)

Adapters/Fittings/Dynacon Style



#69-20310 Dynacon nitrogen coupler x 1/4" hose barb



#69-20311 Dynacon nitrogen coupler x 1/4 - 18 male
#69-21239 Dynacon nitrogen coupler x D.I.S.S. female nitrogen (not shown)



#32-11-20-0011 Replacement knob for Dynacon nitrogen outlets and nitrogen outlets and nitrogen couplers



#69-20584 Dynacon nitrogen adapter x 1/4 - 18 male



#69-20427 Dynacon nitrogen adapter x 1/4" hose barb



#69-20518 Converter, converts Schrader male to Dynacon male - nitrogen (w/set screw lock) - Converter works on swivel adapters only.

Adapters/Fittings

The Industry's only D.I.S.S. Quick Connect

The patented Zip Nut® is the industry's only available D.I.S.S. quick connect. The segmented thread sections allow the nut to be pushed on. Then a 1/4 turn locks & seals the connection. Zip Nut® can be used anywhere a standard D.I.S.S. nut & nipple combination is used.







1) Push on threads

2) 1/4 turn seals it

Select inlet connection required from the left hand column

Select gas service requirements from top row



Inlet	Oxygen	Air	Vacuum	N ₂ O	WAGD	Nitrogen
 1/8" MNPT	69-21398	69-21397	69-21561	69-21562	N/A	69-21563
 1/4" Hose *5/16" Hose	69-21385	69-21383	69-21384	69-21381	69-21494	69-21382
 D.I.S.S. Male	69-21395	69-21393	69-21394	69-21391	69-21396	69-21392
 1/4" MNPT	N/A	69-21565	69-21566	69-21567	69-21568	69-21569

Adapters/Fittings/High Pressure



#69-20451 Oxygen swivel yoke assy x 1/4 -18 male
 #69-20715 Medical air swivel yoke assy x 1/4 - 18 male




#69-20723 Nitrous oxide/oxygen yoke assy x 1/4 - 18 male (N₂O 47.5 - 52.5%)
 #69-20724 Medical diagnostic gas mixture yoke assy x 1/4 -18 male
 #69-20333 Oxygen yoke assy x 1/4 - 18 male
 #69-20334 Nitrous Oxide yoke assy x 1/4 - 18 male
 #69-20335 Medical air Yoke assy x 1/4 - 18 male
 #69-20336 Nitrogen yoke assy x 1/4 - 18 male
 #69-20337 Carbon Dioxide yoke assy x 1/4 - 18 male
 #69-20717 Helium/Oxygen yoke assy x 1/4 - 18 male (HE>80%)
 #69-20718 Helium/Oxygen yoke assy x 1/4 - 18 male (HE<80%)
 #69-20719 Carbon Dioxide/Oxygen yoke assy x 1/4 - 8 male (CO₂<7%)

**Existing large cylinder regulator to be used with D/E cylinders.
 Allows existing large cylinder regulator to be used on small D/E cylinders**



#69-20361 Oxygen Yoke Assy to large cylinder connection
 #69-20452 Medical air yoke assy to large cylinder connection
 #69-20453 Nitrous Oxide yoke assy to large cylinder connection
 #69-20454 Nitrogen Yoke assy to large cylinder connection.



#69-20456 Oxygen yoke assy x 1/4 - 18 male with check valve
 #69-20457 Nitrous oxide yoke assy x 1/4 - 18 male with check valve



#69-20339




#69-21493



#69-20338



#69-21380



#69-20720



#69-21258

#69-21380 Replacement tee handle with plastic grip wheel
 #69-21412 Replacement tee handle with 6" grip (not shown)
 #69-21258 Replacement tee handle with built in "key" accepting cylinder wrench.
 Eases tightening. (1 3/4" long, 16 threads per inch)
 #69-21493 Replacement rubber yoke washer 10/pkg
 #69-20338 Replacement tee handle for standard yokes (1 3/4" long, 16 threads per inch)
 #69-20461 Replacement tee handle (2 3/16" long, 16 threads per inch)
 #69-20462 Replacement tee handle (1 3/4" long, 24 threads per inch)
 #69-20339 Replacement plastic yoke washers 25/pkg (single use only)
 #69-20720 Replacement yoke adapter x 1/4 - 18 male
 #69-20721 Replacement yoke adapter with check valve x 1/4 -18 male (not shown)
 #69-20722 Replacement yoke adapter x 1/2" NPT male (not shown)

Adapters/Fittings/High Pressure

FLEXIBLE STAINLESS STEEL BRAID PIGTAILS*

(3,000 PSI WORKING PRESSURE; 12,000 PSI BURST PRESSURE)

PART NO. DESCRIPTION




#69-20477	Nitrous oxide Large cylinder Pigtail - 18"
#69-20478	Nitrous oxide Large cylinder Pigtail - 24"
#69-20479	Nitrous oxide Large cylinder Pigtail - 36"
#69-20480	Nitrous oxide Large cylinder Pigtail - 48"
#69-20481	Medical air large cylinder Pigtail - 18"
#69-20482	Medical air large cylinder Pigtail - 24"
#69-20483	Medical air large cylinder Pigtail - 36"
#69-20484	Medical air large cylinder Pigtail - 48"
#69-20485	Carbon dioxide Large cylinder Pigtail - 18"
#69-20486	Carbon dioxide Large cylinder Pigtail - 24"
#69-20487	Carbon dioxide Large cylinder Pigtail - 36"
#69-20488	Carbon dioxide Large cylinder Pigtail - 48"
#69-20489	Nitrogen large cylinder Pigtail - 18"
#69-20490	Nitrogen large cylinder Pigtail - 24"
#69-20491	Nitrogen large cylinder Pigtail - 36"
#69-20492	Nitrogen large cylinder Pigtail - 48"
#69-20493	Pigtail with 1/4 - 18 female ends - 18"
#69-20494	Pigtail with 1/4 - 18 female ends - 24"
#69-20495	Pigtail with 1/4 - 18 female ends - 36"
#69-20496	Pigtail with 1/4 - 18 female ends - 48"

See following page for Insufflator fittings.

* Pigtails are also available with check valves. To order check valves, add CV as suffix to product number desired, to receive on inlet side. For check valve on outlet side, add RCV as suffix to product number desired. Flexible pigtails not recommended for oxygen per NFPA99.

Rigid Chrome Plated Copper Pigtails (18" to 20")

Other configurations available with check valves or shut off valves.



#69-20497	Oxygen pigtail
#69-20498	Nitrous oxide pigtail
#69-20499	Medical air pigtail
#69-20752	Oxygen pigtail with check valve
#69-21467	Carbon dioxide pigtail
#69-21468	Nitrogen/Helium tee

Chrome Plated Brass Coupler Tees



#69-20504	Oxygen tee
#69-20505	Nitrous oxide tee
#69-20506	Medical air tee
#69-20507	Carbon dioxide tee (not chrome plated)
#69-20508	Nitrogen/Helium tee
#69-20753	Oxygen tee with check valve




Cap And Chain Assembly

#69-20509	Oxygen female cap and chain
#69-20510	Nitrous oxide female cap and chain
#69-20511	Medical air female cap and chain

Adapters/Fittings/High Pressure

USED WITH INSUFFLATOR SYSTEMS



PART NO.

DESCRIPTION

#69-21202

Carbon dioxide Insufflator yoke, 1/4" SAE Flare with 7/16" - 20 thread




#69-21209

Insufflator pigtail, 36" length 1/4" SAE Flare with 7/16" - 20 thread, both ends




#69-21269

Adapter, 1/4 SAE Flare Female 7/16" - 20 thread x 1/4" hose barb




#69-21267

Adapter, 1/4 SAE Flare Male 7/16" - 20 thread x 1/4" NPT Female



#69-21268

Adapter, 1/4 SAE Flare Male 7/16" - 20 thread x 1/4" NPT Male



#69-21270

Adapter, 1/4 SAE Flare Female 7/16" - 20 thread x 1/4" NPT Male

Adapters/Fittings/High Pressure



PART NO.	DESCRIPTION
#69-21478	Medical mixtures tank nipple with filter
#69-20463	Oxygen tank nipple with filter
#69-20464	Nitrous oxide tank nipple with filter
#69-20465	Medical air tank nipple with filter
#69-20466	Carbon dioxide tank nipple with filter
#69-20467	Helium/nitrogen tank nipple with filter
#69-20468	Oxygen tank nut (CGA 540)
#69-20469	Nitrous oxide tank nut (CGA 326)
#69-20470	Medical air tank nut (CGA 346)
#69-20471	Carbon dioxide tank nut (CGA 320)
#69-20472	Helium/nitrogen tank nut (CGA 580)
#69-21479	Medical mixtures tank nut (CGA 500)
#69-21560	Replacement washer for CO ₂ tank nipple

LARGE CYLINDER CONNECTIONS



BRASS- NOT CHROME PLATED

#69-20731	Oxygen male valve outlet x 1/4 -18 female
#69-20732	Nitrous oxide male valve outlet x 1/4 - 18 female
#69-20733	Medical air male valve outlet x 1/4 - 18 female
#69-20734	Carbon dioxide male valve outlet x 1/4 - 18 female



BRASS- NOT CHROME PLATED

#69-20736	Oxygen male valve outlet outlet x 1/4 - 18 male
#69-20738	Medical air male valve outlet x 1/4 - 18 male
#69-20739	Carbon dioxide male valve outlet x 1/4- 18 male
#69-20740	Helium/nitrogen male valve outlet x 1/4 - 18 male
#69-20737	Nitrous oxide male valve outlet x 1/4 - 18 male



#69-21480	Oxygen tank handnut
-----------	---------------------



#69-21542	One-way check valve, 1/4 - 18 NPT male outlet. Supplied with nitrile o-ring (not for oxygen service)
#69-21543	One-way check valve, 1/4 - 18 NPT female inlet x 1/4 - 18 NPT male outlet Supplied with viton o-ring, for use with oxygen service

Adapters/Fittings

PART NO. DESCRIPTION



- #32-11-30-0004 Lock out adapter for oxygen D.I.S.S. outlets
- #32-11-30-0005 Lock out adapter for vacuum D.I.S.S. outlets
- #32-11-01-0035 Lock out adapter for Chemetron oxygen quick connect outlets
- #32-11-01-0036 Lock out adapter for Chemetron vacuum quick connect outlets
- #32-11-01-0041 Lock out adapter for Chemetron medical air quick connect outlets



- #69-21469 Oxygen extension manifold. Two D.I.S.S. oxygen check valves with wall mount flanges. Hose (2') terminates in 1/8-27 female elbow
- #69-21470 Oxygen extension manifold. Three D.I.S.S. oxygen check valves with wall mount flanges. Hose (2') terminates in 1/8-27 female elbow
- #69-21471 Oxygen extension manifold. Four D.I.S.S. oxygen check valves with wall mount flanges. Hose (2') terminates in 1/8-27 female elbow
- #69-21472 Medical air extension manifold. Two D.I.S.S. medical air check valves with wall mount flanges. Hose (2') terminates in 1/8-27 female elbow
- #69-21473 Medical air extension manifold. Three D.I.S.S. medical air check valves with wall mount flanges. Hose (2') terminates in 1/8-27 female elbow
- #69-21474 Medical air extension manifold. Four D.I.S.S. medical air check valves with wall mount flanges. Hose (2') terminates in 1/8-27 female elbow.
- #69-21475 Vacuum extension manifold. Two D.I.S.S. vacuum check valves with wall mount flanges. Hose (2') terminates in 1/8-27 female elbow.
- #69-21476 Vacuum extension manifold. Three D.I.S.S. vacuum check valves with wall mount flanges. Hose (2') terminates in 1/8-27 female elbow.
- #69-21477 Vacuum extension manifold. Four D.I.S.S. vacuum check valves with wall mount flanges. Hose (2') terminates in 1/8-27 female elbow.













- #69-20346 Oxygen large cylinder test gauge
- #69-20579 Helium or nitrogen large cylinder test gauge
- #69-20580 Carbon Dioxide/oxygen (CO₂ < 7.5%);
Helium/Oxygen (He < 80.5%) CGA-280;
Nitrogen/oxygen (O₂ > 23.5%);
Nitrous oxide/Oxygen (N₂O 47.5% to 52.5%) large cylinder test gauge



- #69-20347 Oxygen in-line pressure test gauge, D.I.S.S. connections
- #69-20348 Nitrous oxide in-line pressure test gauge, D.I.S.S. connections
- #69-20349 Medical air in-line test gauge, D.I.S.S. connections
- #69-20350 Vacuum in-line pressure test gauge, D.I.S.S. connections
- #69-20582 Nitrogen in-line pressure test gauge, D.I.S.S. connections

Adapters/Fittings

	PART NO.	DESCRIPTION
	#69-20650	Plug, 1/8 - 27 male - countersunk
	#69-20651	Plug, 1/4 - 18 male - countersunk
	#69-20957	Cap, 1/8 - 27 female
	#69-20331	Plug, 1/8 - 27 male
	#69-20332	Plug, 1/4 - 18 male
	#69-20540	Close nipple, 1/8 - 27 male x 1/8 - 27 male
	#69-20541	Nipple, 1/8 - 27 male x 1/8 - 27 male - 1 1/2" long
	#69-20542	Nipple, 1/8 - 27 male x 1/8 - 27 male - 2" long
	#69-20543	Nipple, 1/8 - 27 male x 1/8 - 27 male - 2 1/2" long
	#69-20544	Nipple, 1/8 - 27 male x 1/8 - 27 male - 3 1/2" long
	#69-20545	Nipple, 1/4 - 18 male x 1/4 - 18 male - 1 1/2" long
	#69-20546	Nipple, 1/4 - 18 male x 1/4 - 18 male - 2" long
	#69-20547	Nipple, 1/4 - 18 male x 1/4 - 18 male - 2 1/2" long
	#69-20548	Nipple, 1/4 - 18 male x 1/4 - 18 male - 3" long
	#69-20549	Nipple, 1/4 - 18 male x 1/4 - 18 male - 3 1/2" long
	#32-2295-0100	Extension adapter, 1/8 - 23 female x 1/8 - 27 male - 1" long
	#69-20550	Extension adapter, 1/8 - 27 female x 1/8 - 27 male - 2 1/2" long
	#69-20551	Extension adapter, 1/4 - 18 female x 1/4 - 18 male - 3" long
	#69-20684	Extension block - two (2) 1/8 - 27 female outlets x one (1) 1/8 - 27 female centered rear inlet (front view shown), 5" long
	#69-20685	Extension block - two (2) 1/8 - 27 female outlets x one (1) 1/8 - 27 female offset rear inlet (rear view shown), 5" long
		NOTE: 2 outlets spaced 4-1/2" on center measurement
	#69-20687	"L" block, 1/4 - 18 female x 1/8 - 27 female side port
	#69-20691	"L" block, 3 way - 1/4 - 18 female ports, 3" long
	#69-20697	Adapter, elbow - 1/8 - 27 male x 1/8 - 27 female end and 1/8 - 27 female center port
	#69-20726	Adapter, elbow - 1/8 - 27 female end x 1/4 - 27 male with - 1/8 - 27 female center port

Adapters/Fittings


PART NO. DESCRIPTION



#69-21541 Elbow adapter, 90° 1/8-27 female x 5/16" hose barb

#69-21492 Elbow adapter, 90° 1/8-27 female x 1/4" hose barb

It pivots!



#69-21539 Pivot adapter, 90° 1/8 - 27 female x 1/4" hose barb

#69-21540 Pivot adapter, 90° 1/8 - 27 female x 5/16" hose barb



#69-20176 Coupling - 1/4 - 18 male x 1/4 - 18 male

#69-20202 Coupling - 1/8 - 27 male x 1/8 - 27 male


#69-20203 Coupling - 1/8 - 27 male x 1/4 - 18 male



#69-20205 Coupling - 1/8 - 27 female x 1/8 - 27 female

#69-20206 Coupling - 1/8 - 27 female x 1/4 - 18 female

#69-20207 Coupling - 1/4 - 18 female x 1/4 - 18 female



#69-20179 Adapter, 1/8 - 27 female x 1/4" hose barb

#69-20180 Adapter, 1/4 - 18 female x 1/4" hose barb

#69-21523 Adapter, 1/4 - 18 female x 5/16" hose barb



#69-21522 Adapter, 1/8 - 27 male x 5/16" hose barb

#69-20199 Adapter, 1/8 - 27 male x 1/4" hose barb

#69-20200 Adapter, 1/4 - 18 male x 1/4" hose barb

#69-21521 Adapter, 1/4 - 18 male x 5/16" hose barb



#69-20444 Bubble hose nipple, 1/4" hose barb x 1/8 - 27 male



#69-20681 Hose splicer 1/8" hose barb x 1/8" hose barb


#69-20249 Hose splicer 3/16" hose barb x 3/16" hose barb

#69-20250 Hose splicer 1/4" hose barb x 1/4" hose barb



#69-20208 Adapter, 1/8 - 27 female x 1/4 - 18 male

#69-20209 Adapter, 1/4 - 18 female x 1/8 - 27 male
















#69-21201 Bulk head assembly x 1/8 - 27 female



#69-20986 Drop ear, tee 1/8 - 27 female x 1/8 - 27 female x 1/8-27 female

#69-20987 Drop ear, elbow 1/8 - 27 female x 1/8 - 27 female

Adapters/Fittings/Miscellaneous Elbows, Tees, Wyes

	PART NO.	DESCRIPTION
	#69-20682	Adapter - elbow, 1/8 - 27 male x 1/4" hose barb
	#69-20683	Adapter - elbow, 1/4 - 18 male x 1/4" hose barb
	#69-20212	Elbow, 1/8 - 27 female x 1/8 - 27 male
	#69-20213	Elbow, 1/4 - 18 female x 1/8 - 27 male
	#69-20214	Elbow, 1/4 - 18 female x 1/4 - 18 male
	#69-20702	Elbow, 1/8 - 27 female x 1/8 - 27 male, 45°
	#69-20703	Elbow, 1/4 - 18 female x 1/4 - 18 male, 45°
	#69-20215	Elbow, 1/8 - 27 female x 1/8 - 27 female
	#69-20216	Elbow, 1/4 - 18 female x 1/4 - 18 female
	#69-20225	Tee, three-way - 1/8 - 27 x 1/8 - 27 x 1/8 - 27 female
	#69-20226	Tee, three-way - 1/4 - 18 x 1/4 - 18 x 1/4 - 18 female
	#69-20445	Branch tee, 1/8 - 27 female x 1/8 - 27 male x 1/8 - 27 female
	#69-20446	Branch tee, 1/4 - 18 female x 1/4 - 18 male x 1/4 - 18 female
	#69-20447	Street tee, 1/8 - 27 female x 1/8 - 27 female x 1/8 - 27 male
	#69-20448	Street tee, 1/4 - 18 female x 1/4 - 18 female x 1/4 - 18 male
	#69-20449	Cross, 1/8 - 27 female
	#69-20450	Cross, 1/4 - 18 female
	#69-20652	WYE BLOCK, female 1/8 - 27 inlet x two female 1/4 - 18 outlet
	#69-20653	WYE BLOCK, female 1/8 - 27 inlet x two female 1/8 - 27 outlet
	#69-20360	WYE BLOCK, female 1/4 - 18 inlet x two female 1/4 - 18 outlet
	#69-20227	WYE, male 1/4 - 18 inlet x two female 1/4 - 18 outlet
	#69-20228	WYE, female 1/8 - 27 inlet x two female 1/4 - 18 outlet
	#69-20229	WYE, D.I.S.S. female oxygen inlet x two female 1/4 - 18 outlet
	#69-20654	WYE, D.I.S.S. female oxygen inlet x two D.I.S.S. male oxygen outlet
	#69-21234	WYE, 1/8-27 NPT female x Double 1/8-27 NPT female
	#69-21235	WYE, 1/8-27 NPT female x Double 1/4-18 NPT female
	#69-21386	Manifold block, with (6) 1/4 - 18 NPT female ports
	#69-21570	Quad block 1/4 - 18 female ports

Adapters/Fittings/Miscellaneous



#69-21220



#69-21221




#69-21249

PART NO.
#69-21249
#69-21220
#69-21221

DESCRIPTION

OHMEDA male vacuum connector x 1/8 - 27 male
OHMEDA trap bottle connector straight x 1/8 - 27 male
OHMEDA trap bottle connector 45° angle x 1/8 - 27 male



#69-21229
#69-21232
#69-21231
#69-21230
#69-21228

Vacuum Black Type on White
Gas Evacuation Black Type on Purple
N₂O Black Type on Blue
Medical Air Black Type on Yellow
O₂ Black Type on Green adhesive backed aluminum plate - 3" w x 1-1/2" h



#69-21257 Pole clamp, 1/8 - 27 female port on side and 1/8 - 27 female port on bottom for poles up to 1" diameter.

#69-21532 Pole clamp, 1/8 - 27 female port on side and 1/8 - 27 female port on bottom for poles up to 1-1/2" diameter.



Adapters/Fittings/Neb Drive & Rescue Flow

NEB DRIVE

Economical convenient adapter that delivers a preset flow from standard wall outlets. A dependable way to administer aerosol treatments without the need of a flowmeter.

Ohmeda (type)



Flowrate	O ₂	Medical Air
6 to 7 LPM	69-20764	69-20782
7 to 8 LPM	69-20765	69-20783
8 to 9 LPM	69-20766	69-20784
10 to 11 LPM	69-20964	69-20965

RESCUE FLOW STOCK #
14-16 LPM 69-20892
24-26 LPM 69-20898

(see Ohmeda style below)

DISS Female



Flowrate	O ₂	Medical Air
6 to 7 LPM	69-20770	69-20788
7 to 8 LPM	69-20771	69-20789
8 to 9 LPM	69-20772	69-20790
10 to 11 LPM	69-20968	69-20969

RESCUE FLOW STOCK #
14-16 LPM 69-20891
24-26 LPM 69-20897

(see Ohmeda style below)

Schrader (type)



Flowrate	O ₂	Medical Air
6 to 7 LPM	69-20776	69-20794
7 to 8 LPM	69-20777	69-20795
8 to 9 LPM	69-20778	69-20796
10 to 11 LPM	69-20972	69-20973

RESCUE FLOW STOCK #
14-16 LPM 69-20896
24-26 LPM 69-20902

(see Ohmeda style below)

Rescue Flow (Ohmeda Style Shown)



RESCUE FLOW

The "Rescue Flow" is designed to save you time and effort when hooking up your resuscitation bag. The "Rescue Flow" gives a restricted flow of 15 LPM or 25 LPM * from a 50 psi oxygen source and is so small it fits easily on your crash cart.

Chemetron (NCG) (type)



Flowrate	O ₂	Medical Air
6 to 7 LPM	69-20764	69-20782
7 to 8 LPM	69-20765	69-20783
8 to 9 LPM	69-20766	69-20784
10 to 11 LPM	69-20984	69-20985

RESCUE FLOW STOCK #
14-16 LPM 69-20892
24-26 LPM 69-20898

(see Ohmeda style below)

Puritan-Bennett (type)



Flowrate	O ₂	Medical Air
6 to 7 LPM	69-20764	69-20782
7 to 8 LPM	69-20765	69-20783
8 to 9 LPM	69-20766	69-20784
10 to 11 LPM	69-20984	69-20985

RESCUE FLOW STOCK #
14-16 LPM 69-20892
24-26 LPM 69-20898

(see Ohmeda style below)

Oxequip (type)



Flowrate	O ₂	Medical Air
6 to 7 LPM	69-20764	69-20782
7 to 8 LPM	69-20765	69-20783
8 to 9 LPM	69-20766	69-20784
10 to 11 LPM	69-20984	69-20985

RESCUE FLOW STOCK #
14-16 LPM 69-20892
24-26 LPM 69-20898

(see Ohmeda style below)

Med*Star (type)



Flowrate	O ₂	Medical Air
6 to 7 LPM	69-20976	69-20980
7 to 8 LPM	69-20977	69-20981
8 to 9 LPM	69-20978	69-20982
10 to 11 LPM	69-20979	69-20983

RESCUE FLOW STOCK #
14-16 LPM 69-20984
24-26 LPM 69-20985

(see Ohmeda style to left)

*High flow adapter (25 LPM) recommended for use with Laerdal resuscitators & Mercury Disposable CPR Bags.

Adapters/Fittings

VALVES



#69-21222 Wye block w/shut off
1/8 - 27 female inlet x (2) 1/8 - 27 female outlets



#69-21223 1/4 - 18 female x 1/4 - 18 female
#69-21224 1/8 - 27 female x 1/8 - 27 female



NPT Male to Male

#69-21203 1/8 - 27 male inlet, 1/8 - 27 male outlet
#69-21204 1/4 - 18 male inlet, 1/4 - 1/8 male outlet



NPT Female to Female

#69-21205 1/8 - 27 female inlet, 1/8 - 27 female outlet
#69-21206 1/4 - 18 female inlet, 1/4 - 18 female outlet



NPT Male to NPT Female

#69-21207 1/8 - 27 male inlet, 1/8 - 27 female outlet
#69-21208 1/4 - 18 male inlet, 1/4 - 18 female outlet

Compact Valves
250 psi max.

Featuring 1/4 turn shutoff for quick and visual on/off status. Full flow construction for optimum performance.
Durable, chrome plated-brass body. Non-corrosive service.

VALVES

BALL SEAT SHUT - OFF VALVES

Chrome plated brass bodies, Teflon packing, and stainless steel stems with hardened steel balls for efficient, durable use. Large plastic knurled knob for easy operation. *For pressures up to 200 PSI. **For pressures up to 3,000 PSI.



- #69-20623* Valve, 1/4 - 18 male inlet x D.I.S.S. male oxygen outlet
- #69-20624** Valve, 1/8 - 27 male inlet x 1/8 - 27 male outlet
- #69-20625** Valve, 1/8 - 27 female inlet x 1/8 - 27 female outlet
- #69-20626** Valve, 1/4 - 18 male inlet x 1/4 - 18 male outlet
- #69-20627** Valve, 1/4 - 18 male inlet x 1/4 - 18 female outlet



- #69-20680* Dual valve block, D.I.S.S. oxygen female x (2) D.I.S.S. oxygen male outlet w/ valves



- #69-20256* Valve, D.I.S.S. oxygen female inlet x D.I.S.S. male oxygen outlet
- #69-20257** Valve, 1/4 - 18 male inlet x 1/4 - 18 male outlet
- #69-20258** Valve, 1/8 - 27 male inlet x 1/8 - 27 male outlet
- #69-20259** Valve, 1/4 - 18 male inlet x 1/4 - 18 female outlet
- #69-20260** Valve, 1/4 - 18 female inlet x 1/4 - 18 male outlet
- #69-20261** Valve, 1/8 - 27 female inlet x 1/8 - 27 female outlet
- #69-20655* Valve, D.I.S.S. oxygen male inlet x 1/4 - 18 male outlet
- #69-20656* Valve, 1/4 - 18 male inlet x D.I.S.S. oxygen male outlet
- #69-20657* Valve, D.I.S.S. oxygen male inlet x D.I.S.S. oxygen male outlet




- #69-20400 Quick shut-off valve for vacuum, 1/8 - 27 female x 1/8 - 27 male

Stainless Steel Cartridge Valves


(REPLACEABLE CARTRIDGES)

- 316 Stainless Steel Body & Cartridge Seals
- Teflon® TFE & 316 SS Packing
- Machined KEL -F® Soft Tip
- Teflon® TFE Cartridge Seals
- -65° to 250° F



Choose From Three Seat Styles!

Each offering positive shut-off. Seat style K delivers soft-off action with minimal force. Seat style M provides a fine stem movement for accurate gas flow metering. Seat style S is recommended for firm shut-off in corrosive gas applications.



Styles M and S features T-Bar Handle

1/4 - 18 male x 1/4 - 18 male

- #69-20671 Soft seat type K, 3000 PSI
- #69-20672 Metering seat Style M, 6000 PSI
- #69-20673 Metal-to-metal seat style S, 6000 PSI

1/4 - 18 male 90° x 1/4 - 18 male

- #69-20674 Soft seat type K, 3000 PSI
- #69-20675 Metering seat style M, 6000 PSI
- #69-20676 Metal-to-metal seat style S, 6000 PSI

1/4 - 18 female x 1/4 - 18 female

- #69-20677 Soft seat type K, 3000 PSI
- #69-20678 Metering seat style M, 6000 PSI
- #69-20679 Metal-to-metal seat style S, 6000 PSI

Adapters/Fittings/Wye Adapters

MISCELLANEOUS WYE CONFIGURATIONS



- #69-20903** Oxygen duplex adapter - Chemetron male oxygen inlet x two (2) D.I.S.S. male oxygen check valves.
- #69-20904** Medical air duplex adapter - Chemetron male medical air inlet x two (2) D.I.S.S. male medical air check valves.
- #69-20905** Vacuum duplex adapter - Chemetron male vacuum inlet x two (2) D.I.S.S. male vacuum check valves.
- #69-20906** Oxygen duplex adapter - OHMEDA male oxygen inlet x two (2) D.I.S.S. male oxygen check valves
- #69-20907** Medical air duplex adapter - OHMEDA male medical air inlet x two (2) D.I.S.S. male medical air check valves
- #69-20908** Vacuum duplex adapter - OHMEDA male vacuum inlet x two (2) D.I.S.S. male vacuum check valves
- #69-20909** Oxygen duplex adapter - Oxequip 07 male oxygen inlet x two (2) D.I.S.S. male oxygen check valves
- #69-20910** Medical air duplex adapter - Oxequip 07 male medical air inlet x two (2) male medical air check valves
- #69-20911** Vacuum duplex adapter - Oxequip 07 male vacuum inlet x two (2) D.I.S.S. male vacuum check valves.
- #69-20912** Oxygen duplex adapter - Oxequip Med*Star male oxygen inlet x two (2) D.I.S.S. male oxygen check valves
- #69-20913** Medical air duplex adapter - Oxequip Med*Star male medical air inlet x two (2) D.I.S.S. male check medical air valves
- #69-20914** Vacuum duplex adapter - Oxequip Med*Star male vacuum inlet x two (2) D.I.S.S. male vacuum check valves
- #69-20915** Oxygen duplex adapter - Puritan-Bennett male oxygen inlet x two (2) D.I.S.S. male oxygen check valves
- #69-20916** Medical air duplex adapter - Puritan-Bennett male medical air inlet x two (2) D.I.S.S. male medical air check valves
- #69-20917** Vacuum duplex adapter - Puritan-Bennett male vacuum inlet x two (2) D.I.S.S. male vacuum check valves
- #69-20918** Oxygen duplex adapter - Schrader male oxygen inlet x two (2) D.I.S.S. male oxygen check valves
- #69-20919** Medical air duplex adapter - Schrader male medical air inlet x two (2) D.I.S.S. male medical air check valves
- #69-20920** Vacuum duplex adapter - Schrader male vacuum inlet x two (2) D.I.S.S. male vacuum check valves