

Floor Pedestals

STAINLESS STEEL STATIONARY CATHETERIZATION LAB SERVICE PEDESTAL

FEATURES:

- Eliminates hoses and wires that clutter floors.
- Provides easy access to medical gas and electrical services near the point of use, where they are needed.
- Strong 18 gauge type 304 Stainless Steel shell.
- Suitable for remodeling as well as new construction.
- Available with up to 10 services in combinations of medical gas, electric, ground and remote monitor receptacles.

SPECIFICATIONS:

Where indicated on the drawings, provide Chemetron Stainless Steel Catheterization Laboratory Pedestal having 18 gauge sides and 14 gauge bottom plate with No. 4 Satin finishes. All welds are to be ground smooth for seamless appearance. An Access Panel shall be furnished in each pedestal to facilitate installation and servicing of internal fixtures.

Pedestals are available in two sizes: 8"x8"x10" with up to 6 services and 12"x12"x10" with up to 10 services. A stainless steel collar shall be provided to trim out the pedestal against the finished floor.

Factory-installed Medical Gas Outlets shall be pre-piped within the pedestal, and extended 6" beyond the floor mounting plate. Final connection to the Medical Gas Distribution System shall be made by the installer. Medical Gas Outlets shall be Chemetron quick-connect or D.I.S.S. as specified under Medical Gas Outlet Stations and as listed in the Table of Accessories.

The Electrical Contractor shall provide electrical wiring to factory-installed electrical receptacles as indicated in the Table of Accessories.

