

Lamps

LARYNGOSCOPES BLADE BULB SIZE CHART							
TYPE - SIZE	0	1	1 1/2	2	3	3 1/2	4
MACINTOSH	-	SM	-	SM	LG	LG	LG
MILLER	SM	SM	-	SM	SM	-	LG
WISCONSIN	-	SM	SM	LG	LG	-	LG
PHILLIPS	-	LG	-	LG	-	-	-

Standard Small

#10-71907 /Each

Laryngoscope lamp, frosted with knurled grip. For all infant size laryngoscope Blades, Bizzarri-Guiffrida small Blades, Gabuya Blades, MacIntosh Blades - Size 1 and 2, Miller Blades - Sizes 0,1,2 and 3, and Searles and Siker Blades



LARYNGOSCOPE LAMPS

Pre-Focused

#10-71901 5/Pkg (Small)

#10-71902 5/Pkg (Large)

Standard Large

#10-71906 /Each

Laryngoscope lamp, large, frosted with knurled grip. For Bennett Blades, Bizzarri-Guiffrida medium and large Blades - Sizes 3 and 4; Miller Blade - Size 4, MacIntosh Blade - Size 2.5 and 4; Polio and Whitehead Blades, Flagg, Guedel and Wis-Foregger Blades - Sizes 2,3, and 4, and Mercury's PHILLIPS LARYNGOSCOPE® Blades, Sizes 1 and 2 and Robert Shaw.





THE ULTRALITE LARYNGOSCOPE LAMP








Upgrade your present conventional laryngoscope to give a brighter light at greater economy than the fiber optic instrument. Produces a white, bright light up to 2 1/2 times the illumination of conventional lamps. Fits most laryngoscope blades. Made in U.S.A.

#DX-6000 Ultralamp, small (6/pkg)

#DX-6000X Ultralamp, large (6/pkg)

FOR PROPPER INSTRUMENTS		
	Propper No.	Mercury® No.
	#199126	#10-71975
	#199036	#10-71976
FOR FOREGGER INSTRUMENTS		
	Foregger No.	Mercury® No.
	#454003	#10-71907
	#454004	#10-71906

















FIBER OPTIC LAMPS		
	Upsher® (Black Insulator)	#10-19080
	Penlon (Round Handle)	#13-0301015 5/pk
Not Shown	Upsher (Blue Insulator)	#10-20085
	Mercury 2001 Fiberlight®	#10-77301
	Rusch Lynchs and Green Spec	#12-008322 6/pk
	RUSCH FOCS	#12-005630 6/pk

Lamps


























MERCURY MEDICAL'S LAMP CROSS REFERENCE CHART

STATHAM (NATIONAL) INSTRUMENTS #15021 #15022 #15033 #15072 #15073 #15406 #15413	MERCURY ITEM NUMBER #10-719-33 #10-719-34 #10-719-35 #10-719-36 #10-719-37 #10-719-40 #10-719-41
PILLING INSTRUMENTS #49-1500	MERCURY ITEM NUMBER #10-719-42
BAUSCH & LOMB INSTRUMENTS #71-71-44	MERCURY ITEM NUMBER #10-719-46
E.S.I. INSTRUMENTS #49-0002 #49-0003 #49-0004 #49-0010 #49-0026	MERCURY ITEM NUMBER #10-719-47 #10-719-48 #10-719-49 #10-719-50 #10-719-51
A.C.M.I. #49L #51 #52 #54 #56 #57 #58 #58L #58S	MERCURY ITEM NUMBER #10-719-55 #10-719-57 #10-719-58 #10-719-59 #10-719-60 #10-719-61 #10-719-62 #10-719-63 #10-719-64
BOEHM INSTRUMENTS #3 #4 #5 #6 #7 #8 #13 #17 #23	MERCURY ITEM NUMBER #10-719-65 #10-719-66 #10-719-67 #10-719-68 #10-719-69 #10-719-70 #10-719-71 #10-719-73 #10-719-74























Lamps/Welch Allyn Type

ITEM NUMBER	LAMPS WELCH ALLYN #	CURRENT INSTRUMENTS	OBSOLETE INSTRUMENTS
#10-71906	<p># 04800</p> 	<p>Standard Laryngoscopes</p>  <p>664 Series Sizes 2-4 690 Series Sizes 3-4 680 Series Sizes 2-4 692 Series Sizes 3-4 634 Series Sizes 2-4 674 Series Sizes 2-4</p>	<p>Standard Laryngoscopes</p>  <p>694 Series Sizes 3-4 684 Series Sizes 2-4 640 Series Size 1 640 Series Sizes 2-4 40415 Transilluminator</p>
#10-71978	<p># 04900</p> 	<p>Ophthalmoscopes</p>  <p>11720 11730 11735</p>	<p>Ophthalmoscopes</p>  <p>11620 11630</p>
#10-71979	<p># 06000</p> 	<p>Laryngoscope Handles</p>  <p>60813 60814 60713 60815</p>	<p>Laryngoscope Handles</p>  <p>60803 60804</p>
#6783-06100	<p># 06100</p> 		<p>Examination Lights</p>  <p>48400 48435</p>
#6783-06200	<p># 06200</p> 	<p>AudioScopes</p>  <p>23300</p>	<p>AudioScopes</p>  <p>23000, 23020, 23040</p>
#6783-06300	<p># 06300</p> 		<p>Examination Lights</p>  <p>44100 44110</p>

Lamps/Welch Allyn Type

ITEM NUMBER	LAMPS WELCH ALLYN #	CURRENT INSTRUMENTS	OBSOLETE INSTRUMENTS
#6783-03900	# 03900 	Ophthalmoscopes  12810	Ophthalmoscopes  11110
#6783-04100	# 04100 		Examination Lights  48400  48410
#6783-04200	# 04200 	Examination Lights  48700	48600 48610 48625 48635
#10-71932	# 04400 	Ophthalmoscopes  11470	Ophthalmoscopes  11400  11411
#6783-04500	# 04500 		Retinoscopes  18010
#6783-04600	# 04600 		Retinoscopes 
#10-71907	# 04700 	Standard Laryngoscopes  664 Series Size 1  690 Series Sizes 1-2  680 Series Sizes 0-1  692 Series Sizes 1-2  634 Series Sizes 0-1  674 Series Size 1	 684 Series Sizes 0-1  694 Series Sizes 1-2

Lamps/Welch Allyn Type

ITEM NUMBER	LAMPS WELCH ALLYN #	CURRENT INSTRUMENTS	OBSOLETE INSTRUMENTS
#6783-01400	# 01400 		PenLite  77800
#10-71917	# 01800 		Retinoscopes  17700
#6783-02000	# 02000 		Headlights  45003
#10-71927	# 02200 		Otosopes  22900
#6783-02300	# 02300 		Ophthalmoscopes  12900
#6783-02500	# 02500 	Headlights  46003	
#6783-02600	# 02600 	Headlights  49003	
#10-71915	# 03000 	Ophthalmoscopes  11710 StrabismoScope™  12400 EpiScope™  47300	Ophthalmoscopes  11600  11605  11610 Retinoscopes  18000