

# Miscellaneous



**#69-21103**  
Teflon Tape. Pipe thread sealant, 576" roll, 1/4" width.



**#69-21110**  
Vacuum Grease. Silicone lubricant for threads, O-rings, gaskets and seals in vacuum systems. 5.3 oz. tube.



**#69-21111**  
Detachol.™ Great for removal of stubborn adhesives. Safe for most surfaces. 4fl. oz.



**#69-21108**  
Thread Paste Compound 3-1/2 oz. tube.




**#69-21104**  
Heat Sink Compound



**#69-21101**  
Liquid Leak Detector. Pinpoints any gas leak as small as .0001cc per second. Safe for use with oxygen. 8 oz. bottle.



**#69-21105**  
Contact Cleaner 12 oz. can




**#69-21106**  
Spray Lubricant 11-12 oz. can (Brand may vary)



**#69-21102\***  
Service Removable Loctite 50cc bottle. (242)  
**#69-21127\***  
Red Loctite (271) 50cc Bottle  
\* NOTE: Not for use with oxygen.



**#69-21107**  
Freeze Spray Component Cooler 15 oz. can



**#BJ-3900**  
Oxygen Thread Repair (DIE)

# Miscellaneous

#09-0042

Lubricant on wick for  
Bird O-Ring Seals  
(Oxygen compatible)



#69-21113  
**Silver Goop**

Formulated for anti-seize on threaded parts of stainless steel and super alloys.

- Temperature to 1500°F (816°C) .
  - Non-melting.
  - Stays in place.
  - Will not drop off red hot surfaces.
  - Will not give off poisonous metal or oxide fumes at red heat.
  - Balanced composition keeps compound between mating surfaces to prevent seizing, regardless of force applied.
  - Lowers take-up torque on threaded parts, nuts, bolts and pipe threads.
  - Resist water washout.
  - Forms bond with surface to which it is applied.
  - Only occasional application needed after parts have been coated the first few times of use.
  - Not recommended for aluminum or magnesium.
- 1 oz. tube

#09-0631

Lubricant for Bird  
O-Ring Seals  
(Oxygen compatible)



Item Not Shown

#09-3884

Vibratite  
30cc bottle used for various mounting screws and electronic trim pots for Bird ventilators. Great for other similar applications.

**FADE-A-DYNETM**

Removes Betadyne®\* and other iodine-based antiseptic stains from hospital environmental surfaces and fabrics. Just spray area and watch stains disappear.  
#0000-WHS101 16 oz. Spray Bottle

\*Registered Trademark of the Purdue Frederick Co.



#69-21116

Krytox. Fluorinated Grease. Recommended by Puritan-Bennett Corp., as the lubricant of choice for Foregger Anesthesia Machines. 2 oz. Tube



#69-21112  
**High Purity Goop**

Anti-seizing compound for use on titanium, stainless steel, steel, aluminum and high-temperature alloys.

- Insert under conditions to 400°F (204°C).
  - Can be used when oily contamination may otherwise be a problem.
  - Non-corrosive to metals.
  - Impervious to moisture, high voltage and thermal cycling.
- 1 oz. tube



#69-21115

Safety Solvent  
Used to remove residual lubricant from parts during cleaning procedures. 12 oz. Can