


# Nerve Stimulators/components



#BO-10023 Microstim

With the implementation of the FDA regulation, CFR 21 Part #895 and #897, certain devices utilizing exposed metal pins (see photo at left) on lead wires, are required to meet this new guideline. Retrofit Kits are available to convert to the new safety type lead wires. Detailed below are various models and the serial numbers affected.

Neuro Technology Microstim (#BO-10023)  
Serial Numbers

- #EA4113625 and below
- #EG4111381 and below
- #EG40001 and below
- #EF40001 and below



#BO-7100 Microstim Plus

Neuro Technology Microstim Plus (#BO-7100)  
Serial Numbers

- #GA4112145 and below
- #GA40001 and below
- #GH40001 and below

### Microstim Retrofit Kits\*

- #BO-412320\* Gray Retrofit Kit includes matching color panel w/safety jacks, safety lead wires and safety ball electrodes.
- #BO-412322 as above, but for Black colored units.
- #BO-412324 as above, but for P.C. Bone (Cream) colored units



DIGISTIM 3

Neuro Technology **DIGISTIM 3 & 3 PLUS Retrofit Kit\***  
#BO-412326 Retrofit Kit

Neuro Technology **DIGISTIM 2 & 2 PLUS Retrofit Kit**  
#BO-412328 Retrofit Kit



DIGISTIM 2

### \* RETROFIT KITS INCLUDE NEW SAFETY LEAD WIRES

\*Kits require unit disassembly and soldering.

### Replacement Safety Lead Wires

- #BO-100057 Safety Lead Wires w/alligator clips
- #BO-100067 Safety Lead Wires w/snaps
- #BO-100262 Safety Ball Electrodes

### LEAD WIRE WITH "BANANNA" TYPE PLUGS (Safety type)

All accessories comply with FDA performance standard 21 CFR, part 898



#19-00-112P

### Safety Adapter Kit

This Safety adapter kit allows the conversion of older models to accept compliant protected lead sets. Does not alter equipment function. Ball electrodes can still be used normally. Comes complete with 2 adapter jacks, allen key for installation and instructions.



#19-00-108P  
50" Leads with snaps



#19-00-107P  
50" Leads with Alligator Clips



#19-00-110P  
50" Extension Leads w/male & Female ends