



Suction Accessories




#32-22-110006

Bottle holder




#32-22-110005

Anodized aluminum multi-level bottle bracket, 3 feet (.914 m) (does not include #32-22-110007).




#32-22-110024

Suction bottle carrier for one-gallon (3.785 l) collection bottles. Carrier has four conductive swivel casters and a low center of gravity.




#32-22-110004

Bottle Holder Bracket



#32-22-110007

Bottle holder w/lock used with #32-22-110005 multi level bottle bracket.



#32-20020150

Canister holder. Holds additional canisters by attaching to side brackets. (12 per carton).



#69-20921

DISS female to DISS male assembly is used between DISS outlet of suction regulator and Chemetron DISS inlet type disposable collection unit to allow replacement of DCU without the need to remove the suction regulator from the wall. Also serves to overcome space problems.



#69-20328

Vacuum DISS Nut & Stem



#69-20922

DISS to Stem assembly. This DISS female to serrated stem assembly allows the user to mount a suction regulator on the "eye" hook of the #32-20-02-0002 floor mount holder. This permits use of the suction regulator on the operating room floor when stem inlet type disposable collection units are preferred.



#32-20-02-0028

Slide Brackett



#32-20-02-0023
Collection Canister Carrier

Carries up to four Handi-Vak II collection canisters of either 1500 ml or 2400 ml size. Durable, chrome-plated steel construction that's easy to clean. No-slag, free-spinning, rotating casters. Comes with a DISS swivel-nut connection for connecting a suction regulator, a 1/4" barb fitting for attaching a vacuum hose, and 2 shut-off valves for regulating the suction to the canisters.



#32-22-11-0002

Suction Bottle Floor Carrier for one gallon collection bottle



#32-22-11-0003

Suction Stand to mount collection bottles and regulators