


Water Traps

Removes most solid liquid particles from compressed air. Filter element is 5 micron. Available in either manual or automatic drain. Supplied with D.I.S.S. medical air fittings except where noted.



Flow Direction

Blender Application Configurations

- #69-50010 Manual drain type
- #69-50001 Automatic drain type

#69-50002
Manual drain type

#69-50003
Automatic drain type



Flow Direction

Blank Traps

#69-50006 Manual Drain Type (1/8-27 female connections)

#69-50007 Automatic Drain Type (1/8-27 female connections)

#69-50008 Manual Drain Type (1/4-18 female connections)

#69-50009 Automatic Drain Type (1/4-18 female connections)



Wall Or Ventilator Applications/Configurations

Ventilator Configurations (for application to inlet of ventilator)

- #69-50011 Manual Drain Type
- #69-50012 Automatic Drain Type



Flow Direction

- #69-50004 Manual Drain Type
- #69-50005 Automatic Drain Type



Flow Direction


Wall Configurations (for application on wall outlet)

- #69-50015 Manual Drain Type
- #69-50016 Automatic Drain Type



Flow Direction

- #69-50017 Manual Drain Type
- #69-50018 Automatic Drain Type



Flow Direction

Repair Parts and Kits

- #69-50100 5 micron filter element
- #69-50101 Automatic drain only
- #69-50102 Manual drain valve
- #69-50103 Bowl w/automatic drain
- #69-50104 Bowl w/manual drain



**SPECIAL CONFIGURATIONS AVAILABLE
WITH QUICK CONNECT
FITTINGS**