

## The First Full Line of CO<sub>2</sub> Detectors

Finally a Choice in CO<sub>2</sub> Detectors for Patients of Any Size



*No Pull Tab*



### Neo-StatCO<sub>2</sub> <Kg®

Item #: 10-55373

For patients within 0.25 kg - 6 kgs body weight

### ADVANTAGES

- The only CO<sub>2</sub> detector available indicated for use on patients below 1 kg (0.25 kg - 6 kgs)
- Detects up to 100 breaths per minute
- 24-hr performance
- Works in 100% humidity
- Vivid breath-to-breath color change
- 3-year shelf life\*

\* from manufacturing date of CO<sub>2</sub> indicator



### Mini StatCO<sub>2</sub>®

Item #: 10-55371

For patients within 1 kg - 15 kgs body weight

### StatCO<sub>2</sub>®

Item #: 10-55370

For patients over 15 kgs body weight

### ADVANTAGES

- Convenient Pull Tab activates the device when ready to use
- Detects up to 50 breaths per minute
- 24-hr performance
- Works in 100% humidity
- Vivid breath-to-breath color change
- 2-year shelf life\* (StatCO<sub>2</sub> Only)
- 3-year shelf life\* (Mini StatCO<sub>2</sub>)

\* from manufacturing date of CO<sub>2</sub> indicator

Blue color indicates no carbon dioxide is present. Yellow color indicates carbon dioxide is present and ET tube is in the proper place.

Local Distributor: