

ONLY ONE

for CPAP & BiLevel Masks

Deluxe Masks w/ Straight Connector

ORDERING INFORMATION

Part #	Description	Packaged
#10-57103	Large adult deluxe contoured mask with straight connector, head harness	5/Box
#10-57104	Medium/Small adult deluxe contoured mask with straight connector, head harness	5/Box
#10-57105	Child deluxe contoured mask with straight connector, head harness	5/Box



Deluxe Masks w/ Straight Connector & Mask Ports

ORDERING INFORMATION

Part #	Description	Packaged
#10-57112	Large adult deluxe contoured mask with straight connector, 2 mask ports, head harness	5/Box
#10-57113	Medium/Small adult deluxe contoured mask with straight connector, 2 mask ports, head harness	5/Box
#10-57114	Child deluxe contoured mask with straight connector, 2 mask ports, head harness	5/Box



Deluxe Masks w/ 90° Elbow, Anti-Asphyxia

ORDERING INFORMATION

Part #	Description	Packaged
#10-57100	Large adult deluxe contoured mask with 90° elbow, anti-asphyxia, head harness	5/Box
#10-57101	Medium/Small adult deluxe contoured mask with 90° elbow, anti-asphyxia, head harness	5/Box
#10-57102	Child deluxe contoured mask with 90° elbow, anti-asphyxia, head harness	5/Box



Deluxe Masks w/ 90° Elbow, Anti-Asphyxia & Mask Ports

ORDERING INFORMATION

Part #	Description	Packaged
#10-57109	Large adult deluxe mask with 90° swivel port, anti-asphyxia, 2 mask ports, head harness	5/Box
#10-57110	Medium/Small adult deluxe mask with 90° swivel port, anti-asphyxia, 2 mask ports, head harness	5/Box
#10-57111	Child deluxe mask with 90° swivel port, anti-asphyxia, 2 mask ports, head harness	5/Box



ONLY ONE

Low FiO₂ Elbow

Low FiO₂ Elbow:

The Flow-Safe® II has an optional Low FiO₂ Elbow that can be used for specific patients that may require low FiO₂ delivery. The elbow is easy to connect and reduces the delivered FiO₂ to approximately 35%.

Intended Use:

The Flow-Safe II Low FiO₂ elbow is intended for use with Flow-Safe II **ONLY**, without nebulizer.



ORDERING INFORMATION

Part #	Description	Packaged
#10-57128	Low FiO ₂ Non-vented Anti-Asphyxia Valve Elbow	10/Box



11300 - 49th Street North • Clearwater, Florida 33762-4807 USA

Tel. (800) 835.6633 • (727) 573.0088 • Fax (800) 990.6375

www.mercurymed.com