

Distributed by:  
**MercuryMedical**<sup>®</sup>

## PowerMATE<sup>®</sup>/PowerMATE<sup>®</sup> 20 amp\*

(Bold bulleted items refer to PowerMATE 20 amp)

The PowerMATE<sup>®</sup> is the only IV pole mountable power tap  
UL approved for use in patient care areas. UL Recognized  
Component to UL 1363A.

### Features:

- Cord management
- Safe and convenient access to additional outlets
- Increased equipment mobility with 15 ft power cord
- Drip guard protection
- Locks to an IV pole or cart leg
- Patent Pending
- 15 ft, hospital grade power cord - 14 gauge, US NEMA 5-15P molded plug (Gray color)
  - **\*US NEMA 5-20P molded plug (Black color)**
- Hospital grade receptacles - US NEMA 5-15R; 2 pole, 3 wire grounding (Gray color)
  - **\*US NEMA 5-20R; 2 pole, 3 wire grounding (Gray color)**
- Circuit Breaker rating - 15 amps, 125 volts
  - **\*20 amps, 125 AC volts**

### Specifications:

<b>Output</b>	120V, 50-60Hz, 12Amp/* <b>120V, 15Amp</b>
<b>Cord Length</b>	15' Hospital Grade/* <b>15' Hospital Grade</b>
<b>Plug</b>	5-15P Hospital Grade
<b>Receptacles</b>	3 Duplex Hospital Grade
<b>Power ON Indicator</b>	Green Lamp Indicator
<b>Mounting Pole Diameter</b>	0.875" to 1.5"/* <b>0.875" to 1.5"</b>
<b>Operating Environment</b>	0 to 50°C, 0 to 95% RH
<b>Storage Environment</b>	-20 to 65°C, 0 to 95% RH
<b>Weight</b>	3.02 lbs. MAX/* <b>4 lbs.</b>
<b>Approximate Dimensions</b>	5.28" wide, 3.27" deep, 5.65" high/ * <b>5.28" wide, 3.27" deep, 5.45" high</b>
<b>Regulatory Approvals</b>	UL recognized component to UL1363A
<b>Warranty</b>	1 Year
<b>Annual Maintenance</b>	<ul style="list-style-type: none"><li>• Inspect for mechanical damage</li><li>• Check lamp indicator and outlet function</li><li>• Inspect safety labels for legibility</li></ul>
<b>Cleaning</b>	Unplug power and surface clean with 70% alcohol solution
<b>Disposal</b>	Dispose per local regulations for non-hazardous material



#10-20025



**20 amp**  
#10-20026

**WARNING:** Connecting equipment that is NOT part of the medical system to the PowerMATE<sup>®</sup> increases the risk of electrical shock.

## IsoMATE<sup>®</sup>

The IsoMATE<sup>®</sup> contains a medical grade isolation transformer to limit patient leakage current to 100 $\mu$ A. Designed for use with electronic or electromedical devices in patient care areas, the IsoMATE<sup>®</sup> is a UL listed device to UL60601-1 and CSA C22.2 No. 601 US Patent Number 7,777,995.

### Features:

- Mounts to an IV pole with a diameter of 0.75" to 2.50"
- 15 ft, hospital grade power cord - 16 gauge, US NEMA 5-15P modeled plug (Gray color)
- System Rating: 120V, 60Hz, 5 Amp
- Dimensions: 8" wide x 9" deep x 8" high
- Weight: 18 lbs
- 2 year warranty

### Specifications:

<b>Output</b>	120V, 60Hz, 5Amp
<b>Cord Length</b>	15' Hospital Grade
<b>Plug</b>	5-15P Hospital Grade
<b>Receptacles</b>	3 Duplex Hospital Grade - 5 Amp MAX
<b>Power ON Indicator</b>	Green Lamp Indicator
<b>Mounting Pole Diameter</b>	0.75" to 2.5"
<b>Operating Environment</b>	0 to 35°C, 0 to 95% RH
<b>Storage Environment</b>	-20 to 65°C, 0 to 95% RH
<b>Weight</b>	18 lbs. MAX
<b>Approximate Dimensions</b>	8" wide, 9" deep, 8" high
<b>Regulatory Approvals</b>	UL60601-1, CSA C22.2 No. 601
<b>Warranty</b>	2 Years
<b>Annual Maintenance</b>	<ul style="list-style-type: none"><li>• Inspect for mechanical damage</li><li>• Check lamp indicator and outlet function</li><li>• Check transformer leakage current</li><li>• Inspect safety labels for legibility</li></ul>
<b>Cleaning</b>	Unplug power and surface clean with 70% alcohol solution

