

Our Mojo® Vented and Non-Vented masks feature a patented AIR°gel cushion, flexible shell with Custom Fit Technology™ and an Adjustable Tension Wheel. Together, they provide an optimal fit with maximum comfort. Fine-tuned fit and function — have it all in a full-face mask.

## Your Mask, Your Way

### AIR°gel

Maximum comfort, minimum weight. Our patented AIR°gel minimizes pressure points and skin irritation.



### Custom Fit Technology

The shapeable shell of each mask can be easily molded to fit individual needs.

**Flex, fit, forget.**

### 360° Rotating Elbow

For optimal tube positioning and easier mobility.



### AIR°gel Forehead Cushion


### Adjustable Tension Wheel

Removes pressure from the bridge of the nose for a custom, comfortable fit.

### Easy-Release Magnetic Headgear Clips

Lightweight, breathable headgear disconnects from mask easily with one hand.

### Vented and Non-Vented Versions Available



Made in the USA

## What makes our Mojo Mask different?

- *Custom Fit Technology* - a unique soft-shell mask with a shapeable wire allowing for optimal flexibility, fit, and comfort.
- The Adjustable Tension Wheel is exclusive to the Mojo. A wide, continuous range of adjustment ensures a custom fit and removes pressure from the bridge of the nose.
- Available in four sizes with magnetic connectors and simple, breathable quick-release headgear.

### About AIR°gel

Patented AIR°gel cushion delivers an effective seal while minimizing pressure points on the skin.



## Custom Fit Technology

BEND IT.



MOLD IT.



SHAPE IT.



### Mojo Masks with Headgear Specifications:

Product Number	Mask Style	Mask Size	Leak Port	Entrainment Valve	Mask Weight	Mask Dimensions
50835	Vented	Small	✓	✓	144 g	3.73"D x 3.78"W x 6.50"H
50836	Vented	Medium	✓	✓	148 g	3.72"D x 3.76"W x 6.75"H
50837	Vented	Large	✓	✓	156 g	3.79"D x 3.76"W x 7.00"H
50838	Vented	Extra Large	✓	✓	161.5g	3.80"D x 4.45"W x 7.00"H
50920	Non-Vented	Small	✗	✗	144 g	3.73"D x 3.78"W x 6.50"H
50921	Non-Vented	Medium	✗	✗	148 g	3.72"D x 3.76"W x 6.75"H
50922	Non-Vented	Large	✗	✗	156 g	3.79"D x 3.76"W x 7.00"H
50923	Non-Vented	Extra Large	✗	✗	161.5g	3.80"D x 4.45"W x 7.00"H



Does not contain natural rubber latex



# Aura®

The Smart Choice for an Off-the-Shelf Custom Fit Nasal Mask.



### Contact us today.

To place your order, call an authorized distributor.

Call 1.800.742.3646  
 Visit [www.sleepnetmasks.com](http://www.sleepnetmasks.com) to learn more.  
 Sleepnet Corporation, 5 Merrill Industrial Drive,  
 Hampton, New Hampshire 03842 USA

