

The Mojo 2 Vented and Non-Vented masks feature a patented AIR^ogel cushion, moldable shell with Custom Fit Technology™ and a Forehead Adjustment Screw. Together, they provide a personalized fit with maximum comfort.

Forehead Adjustment Screw
Removes pressure from the bridge of the nose for personalized fit.

AIR^ogel[®] Forehead Cushion

Vented, Non-Vented, and Non-Vented AAV Versions Available

Quick-Release Magnetic Headgear Clips
Lightweight, breathable headgear disconnects from mask easily with one hand.

Size Indicator
Available in four sizes

S M L XL

AIR^ogel[®]
Maximum comfort, minimum weight. Our patented AIR^ogel minimizes pressure points and skin irritation.



Custom Fit Technology
The shapeable shell of each mask can be easily molded to fit individual needs.
Bend It. Mold It. Shape It.



Oxygen Port
Allows easy delivery of supplemental oxygen.

360° Rotating Elbow
For optimal tube positioning and easier mobility.



Made in the USA of US and imported parts

What makes our Mojo 2 Mask different?

- Custom Fit Technology - a unique soft-shell mask with a shapeable wire allowing for optimal flexibility, fit, and comfort.
- Adjustment Screw - provides a wide, continuous range of adjustment ensuring a custom fit while removing pressure from the forehead.
- Available in four sizes with magnetic connectors and simple, breathable quick-release headgear.

About AIR°gel

Patented AIR°gel cushion delivers an effective seal while minimizing pressure points on the skin.



Custom Fit Technology

BEND IT.



MOLD IT.



SHAPE IT.



Mojo 2 Masks with Headgear Specifications:

Product Number	Style	Size	Leak Port	Entrianment Valve	Mask Weight	Mask Dimensions
50817	Vented	Small	✓	✓	133.8g	3.60"D x 4.30"W x 7.25"H
50818	Vented	Medium	✓	✓	141g	3.60"D x 4.40"W x 7.40"H
50819	Vented	Large	✓	✓	144.5g	3.60"D x 4.50"W x 7.60"H
50820	Vented	Extra Large	✓	✓	146.5 g	3.60"D x 4.55"W x 7.80"H
50912	Non-Vented	Small	✗	✗	133.2g	3.70"D x 4.30"W x 7.25"H
50913	Non-Vented	Medium	✗	✗	140.4g	3.70"D x 4.40"W x 7.40"H
50914	Non-Vented	Large	✗	✗	143.9g	3.70"D x 4.50"W x 7.60"H
50915	Non-Vented	Extra Large	✗	✗	145.9 g	3.70"D x 4.55"W x 7.80"H
50941	AAV Non-Vented	Small	✗	✓	133.7g	3.65"D x 4.30"W x 7.25"H
50942	AAV Non-Vented	Medium	✗	✓	140.9g	3.65"D x 4.40"W x 7.40"H
50943	AAV Non-Vented	Large	✗	✓	144.4g	3.65"D x 4.50"W x 7.60"H
50944	AAV Non-Vented	Extra Large	✗	✓	146.4g	3.65"D x 4.55"W x 7.80"H



Does not contain natural rubber latex

Contact us today.

General inquiries email
info@sleepnetcorp.com or
 call 1 (603) 758-6600

Visit www.sleepnetmasks.com to learn more.

Sleepnet Corporation, 5 Merrill Industrial Drive,
 Hampton, New Hampshire 03842 USA