

ONLY ONE

PowerMATE[®]/PowerMATE[®] 20 amp*

The PowerMATE[®] is the only IV pole mountable power tap UL approved for use in patient care areas. UL Recognized Component to UL 1363A.

Features:

- Cord management
- Safe and convenient access to additional outlets
- Increased equipment mobility with 15 ft power cord
- Drip guard protection
- Locks to an IV pole or cart leg
- Patent Pending
- 15 ft, hospital grade power cord - 14 gauge, US NEMA 5-15P molded plug (Gray color)
 - ***US NEMA 5-20P molded plug (Black color)**
- Hospital grade receptacles - US NEMA 5-15R; 2 pole, 3 wire grounding (Gray color)
 - ***US NEMA 5-20R; 2 pole, 3 wire grounding (Gray color)**
- Circuit Breaker rating - 15 amps, 125 volts
 - ***20 amps, 125 AC volts**

Specifications:

Output	120V, 50-60Hz, 12Amp/* 120V, 15Amp
Cord Length	15' Hospital Grade/* 15' Hospital Grade
Plug	5-15P Hospital Grade
Receptacles	3 Duplex Hospital Grade
Power ON Indicator	Green Lamp Indicator
Mounting Pole Diameter	0.875" to 1.5"/* 0.875" to 1.5"
Operating Environment	0 to 50°C, 0 to 95% RH
Storage Environment	-20 to 65°C, 0 to 95% RH
Weight	3.02 lbs. MAX/* 4 lbs.
Approximate Dimensions	5.28" wide, 3.27" deep, 5.65" high/ *5.28" wide, 3.27" deep, 5.45" high
Regulatory Approvals	UL recognized component to ULI363A
Warranty	1 Year
Annual Maintenance	<ul style="list-style-type: none">• Inspect for mechanical damage• Check lamp indicator and outlet function• Inspect safety labels for legibility
Cleaning	Unplug power and surface clean with 70% alcohol solution
Disposal	Dispose per local regulations for non-hazardous material

WARNING: Connecting equipment that is **NOT** part of the medical system to the PowerMATE[®] increases the risk of electrical shock.



#10-20025



20 amp
#10-20026

*Bold bulleted items refer to PowerMATE 20 amp